



Science **made** smarter

Service manual

# Eclipse 1.1



Item no. 8507843 – D-0116393 | 2025/07



## Table of contents

|   |            |
|---|------------|
| <b>1. Introduction</b>  | <b>3</b>   |
| 1.1 About this manual   | 3          |
| 1.2 About warnings and cautions                               | 3          |
| 1.3 General information                                       | 3          |
| <b>2. Maintenance and disassembly</b>                         | <b>3</b>   |
| 2.1 General maintenance of Eclipse                            | 3          |
| 2.2 How to clean our products                                 | 3          |
| 2.2.1 Cleaning the OAE probe tip                              | 4          |
| 2.3 Disassembling Eclipse                                     | 6          |
| 2.3.1 Opening the cabinet                                     | 6          |
| 2.3.2 Main board  | 7          |
| 2.3.3 Power supply  | 8          |
| <b>3. Calibration</b>   | <b>9</b>   |
| 3.1 Service, adjustment, and repair                           | 9          |
| 3.2 Equipment requirements                                    | 9          |
| 3.2.1 Equipment requirements for acoustic calibration         | 9          |
| 3.3 Calibration of test equipment                             | 10         |
| 3.4 Starting up from the OtoAccess® Database                  | 10         |
| 3.5 Starting up from NOAH                                     | 11         |
| <b>4. Calibration of transducers</b>                          | <b>12</b>  |
| 4.1 OAE calibration   | 13         |
| 4.2 EP15 - EP25 – VEMP-ABRIS calibration                      | 14         |
| 4.2.1 What the stimuli type curves typically look like        | 18         |
| 4.2.2 Installation of research stimuli                        | 19         |
| 4.3 ASSR calibration  | 23         |
| 4.4 Aided Cortical calibration                                | 24         |
| <b>Appendix A - technical specifications</b>                  | <b>i</b>   |
| EP15/EP25/VEMP/Aided  | iii        |
| TEOAE   | vii        |
| DPOAE   | viii       |
| ABRIS   | ix         |
| ASSR  | x          |
| <b>Appendix B - exploded view</b>                             | <b>ix</b>  |
| <b>Appendix C - spare parts &amp; accessories/disposables</b> | <b>x</b>   |
| <b>Appendix D - update news</b>                               | <b>xvi</b> |



# 1. Introduction

## 1.1 About this manual

This manual is for our multifunctional screening and/or diagnostic device, Eclipse. The contents of this manual cover instrument related data, such as disassembly and calibration. Technical specifications, schematics, parts & accessories list along with an exploded view of the device are found as appendices in the back of the manual.

## 1.2 About warnings and cautions

The below warnings, cautions and notices are applied throughout the manual, indicating the level of attention required for a given action:



**WARNING**

The **WARNING** label identifies conditions or practices that may present danger to the patient and/or user.



**CAUTION**

The **CAUTION** label identifies conditions or practices that could result in damage to the equipment.

**NOTICE**

**NOTICE** is used to address practices not related to personal injury.

## 1.3 General information

We continuously strive to improve our products and their performance, hence the specifications in this service manual are subject to change without further notice.

The performance and specifications of our products can only be guaranteed if technical maintenance is conducted routinely every year. Technical maintenance should be carried out by qualified personnel authorized by Interacoustics.

We are happy to receive any inquiries about our products. Our contact details are:

Interacoustics A/S  
Audiometer Allé 1  
5500 Middelfart  
Denmark

Phone: +45 6371 3555  
Mail: [info@interacoustics.com](mailto:info@interacoustics.com)  
Web: [www.interacoustics.com](http://www.interacoustics.com)



## 2. Maintenance and disassembly



WARNING

Make sure that the instrument is disconnected from the power supply before starting this procedure

### 2.1 General maintenance of Eclipse

The performance and safety of the instrument is maintained when the following recommendations for care and maintenance are observed:

1. To maintain electrical safety during the lifetime of the instrument, a safety check must be made regularly according to IEC 60601-1, Class 1, Type BF if used for ABR measurements and IEC60601-1, Class 1, Type B if used for OAE measurements.
2. Before connecting to the mains, be sure that the local mains voltage and the voltage labelled on the instrument correspond.
3. Do not place the instrument next to a heat source of any kind and to ensure proper ventilation, allow sufficient space around the instrument.
4. We recommend at least one annual instrument evaluation. This ensures that the acoustical, electrical, and mechanical properties are correct. The evaluation should be conducted by a Factory trained technician to guarantee proper service and repair.
5. Observe that no damage is present on the mains cable and connector insulation, and that there is no exposure to any kind of mechanical load involving damage.
6. For maximum electrical safety, switch off the mains power when the instrument is not in use.
7. To ensure the reliability of the instrument, periodic biological measurements should be performed on a person with recorded data. This person could be the operator him/herself.
8. If the surface or parts of the instrument are contaminated, use a soft cloth moistened with a mild solution of water and washing-up liquid or similar. Always disconnect the mains power plug during the cleaning process and be careful that no fluid enters the inside of the instrument or the accessories.
9. After each patient examination, proper cleaning ensures that there is no cross contamination via parts between the patients. General precautions must be observed to prevent transmission of disease from one patient to another. Contaminated ear cushions are to be removed from the transducer before cleaning. Frequent cleaning using water may be used, but periodic use of a mild disinfectant may also be used.
10. Great care should be exercised during the handling of earphones and other transducers, as mechanical shock could cause incorrect calibration.

### 2.2 How to clean our products



Before cleaning, always switch off the device and disconnect from the mains power.  
Use a soft cloth, lightly dampened with cleaning solution to clean all exposed surfaces.  
Do not allow liquid to meet the metal parts inside the earphones / headphones.  
Do not autoclave, sterilize, or immerse the instrument or accessories in any fluid.  
Do not use hard or pointed objects to clean any instrument parts or accessories.  
Do not allow parts that have been in contact with fluids dry before cleaning.  
Rubber ear tips or foam ear tips are single-use only components.  
Ensure that isopropyl alcohol does not meet any screens on the instruments.  
Ensure that isopropyl alcohol does not meet any silicone tubes or rubber parts.

**Recommended cleaning and disinfection solutions:**

Warm water with mild, nonabrasive cleaning solution (soap)  
Normal hospital bactericides  
70% isopropyl alcohol

**Procedure:**

Clean the instrument by wiping the outer case with a lint free cloth lightly dampened in cleaning solution.  
Use a lint free cloth dampened in cleaning solution for the cushions, patient tone switch and other parts.  
Make sure to not get moisture in the speaker portion of the earphones and similar parts.



## 2.2.1 Cleaning the OAE probe tip

To record correct OAE measurements, it is important to make sure that the probe system is kept clean at all times.



Do not use pins or threads/needles to remove deeply positioned deposits in the small tubes inside the probe tip, as two of the channels holding acoustic filters may pop out or be damaged. Instead, use the extra probe tip replacement parts supplied together with the OAE system. The probe and cables can be cleaned with alcohol wipes. If the OAE module is used as an infant screener within a hospital setting, the pediatric ward will specify disinfection procedures and recommend the appropriate agents. In this case, the probe should be cleaned after every measurement. Also, consider wiping Eclipse thoroughly.

Never wash or dry the probe tips with a temperature above 70°C/ 158°F.

Do not immerse the probe tip in water.

### NOTICE

Never clean the probe housing by immersion in a solution.

The below illustrates a disassembled OAE probe.



Figure 1

Never attempt cleaning the probe tip while the tip is still fitted to the probe housing. Instead, remove the probe cap by turning it counterclockwise. The probe tip can then be removed from the housing.




Figure 2

Pull the cleaning tool apart to access the thin brush and thin rigid plastic cord inside.




Figure 3

Use the plastic cord or brush to remove debris from the probe tip. Always insert the cleaning tool from larger end of the tip – this prevents accumulation of debris inside the small, round vents. Make sure to insert the cleaning tool into the open channel vent, as the other vents contain small red acoustic filters.



Figure 4

Alternatively, use a Thornton Bridge & Implant Interdental Cleaner to remove debris. Thread the stiff end of the cleaning thread into the rear end of the open probe tip tube and pull the floss through.

When finished, place the probe tip back on the probe. Make sure that the holes fit into the corresponding cavities.

Screw the probe cap back onto the probe. Finger force is sufficient to tighten the screw, hence tools to affix the probe cap is not recommended.



## 2.3 Disassembling Eclipse



WARNING

Make sure that the instrument is disconnected from the power supply before starting this procedure.

The following paragraphs contain the most common actions related to Eclipse, namely, how to disassemble the instrument for service. Carried out by a technician authorized by Interacoustics, the most common action would be replacing the mainboard.

Calibration, another common action that this instrument requires, is described thoroughly in [Chapter 3](#).

Tools required: Torx T10

### 2.3.1 Opening the cabinet

Disassembling Eclipse is a simple procedure of three steps.

Remove the two screws holding the upper cabinet in place.



Figure 5

Then remove the upper cabinet by sliding it towards yourself, away from the rear panel.




Figure 6



Figure 7

#### NOTICE

Before sliding the upper cabinet fully back, the ground cable going from the upper cabinet to the lower cabinet ground connection, must be released.

Otherwise, it hinders the upper cabinet from sliding off.



### 2.3.2 Main board

This chapter describes how the mainboard is removed.




Figure 8

Remove the 5 nuts keeping the mainboard in place.



Figure 9

Remove the nut on the OAE opening.



Figure 10

Remove the 2 screws holding the preamp socket in place.




Figure 11

Finally remove the cables indicated to the left.

The mainboard can now be removed from the lower cabinet.



### 2.3.3 Power supply




Figure 12

Remove the screws holding the power supply in place.  
Remove the cable indicated on the left from its socket.  
The power supply can now be removed from the lower cabinet.



### 3. Calibration

This chapter contains calibration standards for Eclipse. Included here is an explanatory guide and a complete set-up of how to calibrate each transducer type.

#### 3.1 Service, adjustment, and repair

##### Annual calibration

This instrument can be serviced, adjusted, and repaired without losing the CE marking validity, provided that the measuring equipment used is fulfilling the demands below, the adjustment procedures are followed and technicians are authorized by Interacoustics.

To ensure correct stimuli, we recommend a yearly calibration of the following modules:

|                        |  |             |                     |
|------------------------|--|-------------|---------------------|
| TEOAE CALIBRATION      | Probe stimuli are calibrated in peSPL values using the IEC 711 ear simulator coupler made in accordance to IEC 60318-4.                      |             |                     |
| DPOAE CALIBRATION      | Probe stimuli L1 and L2 are calibrated individually in SPL values using the IEC 711 ear simulator coupler made in accordance to IEC 60318-4. |             |                     |
| ABR / ASSR CALIBRATION | IEC 60645-1/ANSI S3.6  | IEC 60645-3 | IEC 60645-7, Type 1 |

#### NOTICE

We recommend that a calibration check label with the next date of calibration, is attached to the equipment.

#### 3.2 Equipment requirements

Calibrating Eclipse encompasses both hardware calibration and transducer calibration. In the following paragraphs, the equipment requirements for both types are listed.

##### 3.2.1 Equipment requirements for acoustic calibration

To ensure proper calibration of audiometric equipment, test instruments must be reliable, stable and their calibration certified. The minimum requirements for measuring equipment are:

- 1) Measuring amplifier with condenser microphone input, or a sound level meter fulfilling the demands of IEC 651 Type 1
- 2) Acoustical coupler in accordance with IEC60318-4 (IEC 711), Occluded ear simulator IEC60318-1 artificial ear and IEC60318-6 artificial mastoid
- 3) An acoustic calibrator
- 4) A general purpose oscilloscope



Figure 13



Figure 14



Figure 15

Appendix B  
contains an  
**exploded  
view**  
of  
Eclipse 1.1



### 3.3 Calibration of test equipment

All test equipment used for the calibration of audiometers is to be calibrated annually. Calibration of all test equipment is to be performed by a laboratory which is traceable to the National Institute of Standards and Technology (NIST) or equivalent. Test equipment calibration documents are to be kept on file.

### 3.4 Starting up from the OtoAccess® Database

Even with no hardware detection, the selected Eclipse module still opens. However, it is not possible to start the actual test before the hardware is properly detected.

Therefore, make sure that Eclipse is switched on and correctly connected before opening the software module.

To get started from the OtoAccess® Database:




Figure 16

For more information, refer to the OtoAccess® Database Instructions for Use



### 3.5 Starting up from NOAH

The following concerns use with ASSR only.

Having opened NOAH, select the current patient's name. For patients not yet listed, click *Add a New Patient*, and fill in the required fields.

Then save the patient details by clicking *OK* and select the ASSR module from the list.




Figure 17

For more information about working with the database, please refer to the operation manual for NOAH



## 4. Calibration of transducers

This chapter contains a transducer calibration with Eclipse, exemplified with an EARTone ABR transducer.

First, set up the calibration equipment as illustrated below. The calibration setup is in fact a circuit.




Figure 18

For the following modules, a yearly calibration check to ensure correct stimuli, is recommended.

|                        |  |             |                     |
|------------------------|--|-------------|---------------------|
| TEOAE CALIBRATION      | Probe stimuli are calibrated in peSPL values using the IEC 711 ear simulator coupler made in accordance with IEC 60318-4.                      |             |                     |
| DPOAE CALIBRATION      | Probe stimuli L1 and L2 are calibrated individually in SPL values using the IEC 711 ear simulator coupler made in accordance with IEC 60318-4. |             |                     |
| ABR / ASSR CALIBRATION | IEC 60645-1/ANSI S3.6  | IEC 60645-3 | IEC 60645-7, Type 1 |



## 4.1 OAE calibration




Figure 19

**TEOAE**

Start the Eclipse OAE Calibration program. Choose *TEOAE* from the main menu, and then *Transducer Calibration*.

Follow the on-screen instructions.




Figure 20

**DPOAE**

Choose *DPOAE* from the main menu, and then *Transducer Calibration*.

Follow the onscreen instructions.

TEOAE calibration is always required, even when no DPOAE is included



## 4.2 EP15 - EP25 – VEMP-ABRIS calibration

Then, switch on Eclipse and open the calibration programme via the below path:




Figure 21





Figure 22


When the calibration program is running, start with the *Stimuli Level* tab, and choose the transducer to calibrate.

First, calibrate with the stimuli type *1KHz Click ref.*

Then observe the values and follow the on-screen instructions.


Repeat for each transducer up for calibration.





This is what the reference tone curve looks like.

Figure 23



Then, repeat the same for the stimuli type *Click Traditional*.

These options are not in use.

Figure 24




Figure 25

This is what the click looks like on the oscilloscope.




Figure 26

Then click the *Tone Bursts, Click, Chirp® & NB Chirps®* tab, and choose the transducer to calibrate.

Use the below table to adjust the stimuli and frequencies correctly.

Follow the on-screen instructions and observe the values.



| Stimuli   | Burst | NB CE-Chirp | Click | CE Chirp |
|-----------|-------|-------------|-------|----------|
| Frequency |       |             |       |          |
| 250 Hz    | X     |             |       |          |
| 500 Hz    | X     | X           |       |          |
| 750 Hz    | X     |             |       |          |
| 1 kHz     | X     | X           | X     | X        |
| 1.5 kHz   | X     |             |       |          |
| 2 kHz     | X     | X           |       |          |
| 3 kHz     | X     |             |       |          |
| 4 kHz     | X     | X           |       |          |
| 6 kHz     | X     |             |       |          |
| 8 kHz     | X     |             |       |          |

This table displays the relevant types of stimuli and their matching frequencies.






**NOTICE**

Calibrating the stimuli CE-Chirp High and CE-Chirp Low is not required.

Figure 27




### 4.2.1 What the stimuli type curves typically look like





### 4.2.2 Installation of research stimuli



Open the *Wave files* tab and click *Add Stimulus...* to import the research stimuli. Also import 1kHz\_warble.wav for calibration.

**NOTICE**

Free field calibration requires a dedicated transducer (e.g. the Left Headset).<sup>1</sup>

Figure 33




Figure 34

<sup>1</sup> When calibrated to free field speakers, the selected output cannot be used for the original transducer. In this example, the Left Headset must be selected in the EP25 software in order to present stimuli to the speaker.




Figure 35

Select 1 kHz warble.wav and adjust to the given level (e.g. 100dB SPL).

Measure the peak-peak voltage on the oscilloscope.

### EQUIPMENT REQUIREMENTS

711 coupler (IEC60318-4) for insert phones

Artificial ear (IEC60318-1) for supra aural phones

Artificial mastoid (IEC60318-6) for bone conductors

Sound level meter

Oscilloscope (connected to the output of the sound level meter)




Figure 36

Then, check the 100 dB peSPL level and adjust the calibration slider until the same Vp-p as the reference tone in the previous steps is reached.





Figure 37

Uncheck the box *Continues playback (when used as masking stimulus)*.



Select 0 dB in the *Level (peSPL)* box.

Then adjust the slider *peSPL<nHL* to 0.5 dB.

Figure 38



### 4.3 ASSR calibration

Open the ASSR calibration software from the file as below:




Figure 39





Figure 40

In the ASSR calibration program, start with the *ASSR Frequencies* tab and choose the transducer to calibrate.

Start with *Ref. Tone* and follow instructions. Follow this guide through *ASSR Frequencies* and *WN*.

Then observe the values and follow the on-screen instructions.

Repeat for each transducer up for calibration.





#### 4.4 Aided Cortical calibration

NOTE: Unlike the other transducers and modules delivered with the Eclipse, the Aided Cortical module and loudspeaker are NOT calibrated from the factory. Hence, a calibration of the Aided Cortical setup is ALWAYS needed before first time use.

For the calibration of the Aided Cortical module, the following calibration equipment is needed.

- A Sound Level Meter (SLM) (i.e. B&K Type 2250 or similar)
- A stand or tripod to position the SLM in the correct position/height

##### Installation:

Before calibration can begin, the equipment needs to be installed. The speaker stand must be assembled and the speaker placed on top (secure it with the supplied finger-screw on the bottom). Adjust the height of the speaker to ensure its center is roughly in the same height as where the patient's head is expected to be. Secure the height with the pin and tighten the knob by hand.



Figure 41



**Setup:**

The speaker is connected via socket “LINE 2” with the supplied cable. The power supply is connected to either of the power sockets. Do not switch on the speaker before it has been connected to the Eclipse.




Figure 42



Figure 43

**Volume:**

It is recommended to turn the overall speaker volume up with the slit positioned vertically (at “12” o’clock) for the initial setup. This will set the volume at roughly 50% and it will be adjusted during the calibration.



Figure 44



**Microphone stand assembly:**

Assemble the microphone stand including the microphone adapter at the top of the pole. The microphone is attached by sliding the cable into the groove and pushing the microphone gently into the holder until fixed.



Figure 45



Figure 46



Figure 47



**Connection:**

To connect the speaker and microphone cable to the Eclipse, the supplied 3.5 mm to 6.5 mm stereo-jack adaptors are used for both cables.

The speaker is connected to the Eclipse via the “Bone” socket and the microphone is connected to the Eclipse via the “Talk Forw.” socket.



Figure 48

The speaker can also be connected using the optional 10 m SP100 cable with one 6.35 mm plug and one 3.5 mm plug.

**Placement:**

Place the speaker in the desired position – it is recommended to use a position away from the walls in the room. Place a chair 1.5 m in front of the speaker. Make sure that there is room for this setup.

During calibration, the patient chair is removed, and the microphone stand is placed in the location where the chair was placed. Adjust the height of the microphone stand to be in the same level as the patient’s ear will be while they are sitting in the chair. Face the tip of the metal tube of the microphone toward the loudspeaker.

Assemble the SLM to the tripod and place the tip of the SLM so it is facing the speaker and parallel with the microphone. The SLM should be in a similar height as the microphone but the two may not touch. The SLM tip must be at the same plane as the tip of the microphone.



Figure 49




Figure 50



Turn on the SLM and set the testing parameters as follows:

- Fast mode
- C frequency weighting – alternatively LIN or Z weighting (do not use A weighting)
- Ideally with large numbers for easy view during the calibration

The free field setup is sensitive to influence from environmental noises like fans, doors slamming, people talking and so forth. Thus, it is a must to minimize or remove any such noise sources during the calibration process.

Calibration software:

Once all parts are placed and turned on (speaker and SLM), open the calibration software on the PC connected to the Eclipse. Choose the tab labeled “Aided Cortical stimuli”.

Note that as soon as you choose this tab, a sound will be played constantly from the loudspeaker!

Once in the tab, the steps are explained in the dialog box.




Figure 51



The first step is to adjust the overall volume of the setup using the volume control of the speaker (on the SP90A, use the volume control on the rear of the speaker to reach the required level).



Figure 52

Note that the “Speaker volume” can be used as a visual guide to ensure the system is working correctly. The displayed volume is the SPL picked up by the microphone near the SLM. The value of the SLM and the value displayed by the “Speaker volume” may not be completely the same. Within 10 dB of each other is acceptable.




Figure 53



Next step is to choose the first stimuli for calibration – the ManU-IRU Low.




Figure 54

Once selected, the stimulus will start playing and the SPL value can be read from the SLM. The target is 80.0 dB. If the value is incorrect, use the slider to either make the stimulus louder or softer. The slider moves in steps of 0.5 dB, so choose the value closest to 80.0 dB.



Once the value is adjusted, click on the ManU-IRU Medium and repeat the process.




Figure 55

Repeat the steps for all the HD-Sounds and LING. Note the "Reserved 1-3" do not contain any stimuli and should be skipped during the calibration.



Once all the stimuli are calibrated, an optional correction factor adjustment is possible. For most applications, these values should NOT be adjusted. Only in special cases where the end-user has specific requirements should these be adjusted. If no adjustment is required, jump to the final step in the calibration process.

The process is similar to the 80 dB calibration – a specific stimulus is selected and the “0 dB” slider is selected. With the slider, the correction factor for the chosen stimulus is adjusted. Once set, the next stimulus is selected and the process is repeated.




Figure 56



The final step in the calibration is to set the reference level for the Sound Field Analyzer. This is done by clicking on “Stimulate and measure” in the bottom section of the calibration tap. Note as soon as this is selected, a Pink Noise signal will be played from the loudspeaker.




Figure 57

The Pink Noise signal will be recorded by the Sound Field Analyzer via the connected microphone and the values for each frequency will stabilize at a fixed level. Small (up to ~1 dB) fluctuations to the dB value are acceptable. Once stable, the calibration is done. Press the “OK” button below to store the results and exit the calibration setup.



There is an option to set the level of the test signal used by the Sound Field Analyzer in the “Settings” tab. It is recommended to use 75 dB. If needed, another level can be chosen.




Figure 58



## Appendix A - technical specifications

|                                 |  |   |
|---------------------------------|--|---|
| <b>Medical CE-mark:</b>         | The CE-mark in combination with MD symbol indicates that Interacoustics A/S meets the requirements of the Medical Device Regulation (EU) 2017/745 Annex I<br>Approval of the quality system is made by TÜV – identification no. 0123 |   |
| <b>Standards:</b>               | <b>Safety:</b>   | IEC 60601-1:2005+AMD1:2012+AMD2:2020 (Class I, Type BF)<br>IEC 60601-2-40:2016<br>CAN/CSA-C22.2 No.60601-1:2014/A2-2022<br>ANSI/AAMI ES60601-1:2005/A2:2021 |
|                                 | <b>EMC:</b>  | IEC 60601-1-2:2014+AMD1:2020  |
| <b>Power Supply:</b>            | <b>Input Volts:</b>  | 100 –240VAC, 50/60Hz.   |
|                                 | <b>Consumption:</b>  | 26W (0.3A Max)  |
|                                 | <b>Safety marking</b>  |    |
| <b>Operating environment:</b>   | <b>Operating Temperature:</b>  | 15 – 35 °C (59 - 95°F)  |
|                                 | <b>Rel. Humidity:</b>  | 30 – 90%  |
|                                 | <b>Ambient Pressure:</b>   | 98kPa – 104kPa  |
| <b>Transport &amp; Storage:</b> | <b>Storage Temperature:</b>  | 0°C – 50°C (32°F - 50°F)  |
|                                 | <b>Transport Temperature:</b>  | -20 – 50 °C (-4°F - 122°F)  |
|                                 | <b>Rel. Humidity:</b>  | 10 – 95% (non condensing)   |
| <b>Warm up time:</b>            |  | 10 minutes at room temperature (20 °C) (68°F).  |
| <b>General</b>                  |  |   |
| <b>PC control:</b>              | <b>USB:</b>  | USB 1.1 or 2.0 for input/output for computer communication. Eclipse if fully operated from a PC   |
| <b>Construction:</b>            |  | Metal cabinet   |
| <b>Eclipse Dimensions</b>       |  | (L x W x H) 28 x 32 x 5.5 cm (11 x 12.6 x 2.2 Inches)   |
| <b>Eclipse Weight</b>           |  | 2.5kg / 5.5 lbs excluding accessories   |



## EP15/EP25/VEMP/Aided

|                                     |  |  |   |
|-------------------------------------|--|--|---|
| <b>Medical CE-mark:</b>             | The CE-mark in combination with MD symbol indicates that Interacoustics A/S meets the requirements of the Medical Device Regulation (EU) 2017/745 Annex I<br>Approval of the quality system is made by TÜV – identification no. 0123 |  |   |
| <b>Standards:</b>                   | <b>Test Signal:</b>  | IEC 60645-3, 2007  |   |
|                                     | <b>AEP</b>   | IEC 60645-7, 2009. Type 1  |   |
| <b>EPA Preamplifier:</b>            | <b>Two channels standard</b>   | EPA4 Cable Collector (4 electrodes). Standard 50 cm. Option: 5 cm or 290 cm  |   |
|                                     | <b>One Channel (optional)</b>  | EPA3 Cable Collector (3 electrodes). 50 cm   |   |
|                                     | <b>Gain:</b>   | 80 dB/60 dB (60dB = VEMP (-20dB))  |   |
|                                     | <b>Frequency response:</b>   | 0.5 – 11.3kHz  |   |
|                                     | <b>CMR Ratio:</b>  | Minimum 100 dB. Typical 120 dB @55 Hz  |   |
|                                     | <b>Noise (RTI)</b>   | =< 15 nV/√Hz   |   |
|                                     | <b>Radio frequency immunity:</b>   | Typically 20 dB improvement over previous available designs  |   |
|                                     | <b>Max input offset voltage:</b>   | 2,5 V  |   |
|                                     | <b>Input impedance:</b>  | >=10 MΩ/=< 170 pF  |   |
|                                     | <b>Power from main unit:</b>   | Insulated power supply with 1500 V isolation. The signal is digitally/capacitive insulated.  |   |
| <b>Specifications as EPA4</b>       | <b>Impedance measurement:</b>  | Selectable for each electrode  |   |
|                                     | <b>Measurement frequency:</b>  | 33 Hz  |   |
|                                     | <b>Waveform:</b>   | Rectangular  |   |
|                                     | <b>Measurement current:</b>  | 19µA   |   |
|                                     | <b>Range:</b>  | 0.5 kΩ – 25 kΩ   |   |
|                                     | <b>Stimulus:</b>   | <b>Stimulus rate:</b>  | 0.1 to 80.1 stimuli per second in steps of 0.1.                 |
| <b>Envelopes/Windows:</b>           |  | Bartlett, Blackman, Gaussian, Hamming, Hanning, Rectangle and Manual (Rise/Fall and Plateau)   |   |
| <b>Masking:</b>                     |  | White noise. Calibrated and presented in peSPL.  |   |
| <b>Transducer:</b>                  |  | Insert phone, calibrated on an IEC 711 coupler.<br>Headphone with independent calibration (optional)<br>Bone conductor (optional)<br>Free field loudspeaker (optional)   |   |
| <b>Level:</b>                       |  | 20 – 135.5 dB peSPL, please refer to chapter stimulus maximum intensity for converted nHL as this depends on the frequency.  |   |
| <b>Polarity:</b>                    |  | Condensation, Rarefaction, Alternating.  |   |
| <b>Click:</b>                       |  | 100 µs (200Hz -11kHz)  |   |
| <b>Tone Burst Frequency:</b>        |  | 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000 and 8000 Hz.   |   |
| <b>Tone Burst Stimulation Time:</b> |  | Stimulation up to 780 ms   |   |
| <b>NB CE-Chirp® LS Freq.:</b>       |  | 500, 1000, 2000 and 4000 Hz  |   |
| <b>Broadband CE-Chirp®: LS</b>      |  | 200Hz -11kHz   |   |
| <b>Relative Masking Level:</b>      |  | +30dB to -40 dB relative to stimulus level. The stimulus level is presented in nHL. The masking level is only presented in SPL, and can therefore not exceed the loudness of the stimulus. E.g. a stimulus presented at 100dBnHL, and relative masking level at 0dB would provide a masking level of 100dB peSPL. This would equal a level 75dBnHL.<br>Maximum masking levels:<br>Insert phones: 110dB SPL, relative levels 0 to -40.<br>Headphones: 110dB SPL, relative levels 0 to -40.<br>Insert phones: 110dB SPL, relative levels +60 to -40. |   |
| <b>Absolute Masking Level:</b>      |  | 0dB to 110 dB SPL absolute level. The masking level is only presented in SPL, and can therefore not exceed the loudness of the stimulus. E.g. a stimulus presented at 100dBnHL, and relative masking level at 0dB would provide a masking level of 100dB peSPL. This would equal a level 75dBnHL.<br>Maximum masking levels:<br>Insert phones: 110dB SPL, relative levels 0 to -40.<br>Headphones: 110dB SPL, relative levels 0 to -40.<br>Insert phones: 110dB SPL, relative levels +60 to -40.   |   |
| <b>Recording:</b>                   |  | <b>Analysis Time:</b>  | -150 ms prior to stimuli and up to 1050 ms (license dependent). |
|                                     |  | <b>A/D Resolution:</b>   | 16 bit.   |
|                                     | <b>Sampling frequency</b>  | 30 kHz   |   |
|                                     | <b>Artifact Reject System:</b>   | Standard voltage based system  |   |
|                                     | <b>Rejection levels:</b>   | Manual 0.2 - 640 µV input with 0.1uV steps.  |   |
|                                     | <b>Anti-aliasing filter:</b>   | Internal filter in ADC   |   |
|                                     | <b>Dots per Trace:</b>   | 450 displayed.   |   |



|                               |                              |   |
|-------------------------------|------------------------------|---|
|                               | <b>Low Pass Filter:</b>      | None or 17 – 12000 Hz, depending on the measurement type. 33 taps FIR Filter without wave peak latency displacement.  |
|                               | <b>High Pass Filter:</b>     | 0.83 Hz to 500 Hz depending on the measurement type.  |
|                               | <b>DSP Low Pass Filter:</b>  | None, 100, 300, 750, 1k, 1.5k, 2k, 3k, 4k, 5k, 7.5k Hz  |
|                               | <b>DSP High Pass Filter:</b> | 0.5, 1.0, 3.3, 10, 33, 100 Hz   |
| <b>Display Gain:</b>          |                              | General Display Gain. Applicable during testing. Single Curve Display Gain. Applicable during testing.  |
| <b>Controlled parameters:</b> |                              | Stimuli Rate, Number of stimuli, Polarity, Click, Tone Burst (Frequency, no. of sine waves, window), Stimulus intensity, Number of curves per intensity, Intensity (Ascending, Descending), Soft attenuator, Stimulus ear, Transducer, Masking level, Preliminary filter setting, Recording onset, Automatic next intensity (Wave repro level on screen), General Display Gain, Single Curve Display Gain, Baseline, Latency norm, Report templates, Print out, Manual stimulus to familiarization, Talk Forward. |
| <b>Data collection:</b>       |                              | Impedance test, Waveform buffer (A/B, Contra, Ipsi-Contra, A-B = Noise), Curve (Hide, Fixate, Merge, Delete), Online EEG, Waveforms storage in unlimited storage database.  |
| <b>Data Recovery:</b>         |                              | Lost data due to crash of Windows® will in almost all cases be available upon re-establishing Windows® operation.   |

**NOTICE**

The transducer settings match the individual Eclipse, as the calibration data is stored on the instrument. Using other transducers require individual calibrations of the transducers.

Two types of calibration are used with AEP:<sup>2</sup>

| peSPL (peak equivalent Sound Pressure Level)  | nHL (normal Hearing Level)   |
|---|--|
| Objective measure of the sound stimulus pressure level.   | Compensating correction for deviation in perceived loudness of the brief stimuli like Clicks and Tone bursts.  |
| For a given peSPL dB value, the max. acoustical or vibration level is calibrated to match the level of continuous tones used to obtain the sound level meter dB SPL reading. As the duration of sound stimuli for AEP are short, the energy delivered is not perceived with the same subjective loudness as the equivalent stimulus would provide, were it a continuous tone. Hence, the acoustical or vibration value given in dB peSPL does not correspond accurately with normal HL figures. For Clicks, there is a 35.5 dB deviation (70 dB peSPL sounds as 35 dB HL), and for Tone Bursts, the deviations are in the 20-30 dB range (depending on frequency and number of sine waves used in burst). The transducers limit the stimulus intensity to 135.5 dB peSPL. | There is a direct similarity between the indicated level in nHL and the HL levels well known from normal audiometry. Brief tone burst correction values from peSPL to nHL are based on a 2-1-2 manual burst as specified in ISO 389-6-2007. Longer duration tone bursts for AMLR and ALR are employing peRETSPL values similar to continuous pure tones (as used in conventional audiometers) since the temporal integration of tones lasting at least 50 ms is considered sufficiently trivial to ignore. |

The calibration values are found on the next page.

<sup>2</sup> Reference calibration values for the CE-Chirp® stimulus family are not specified in the current international calibration standard (ISO 389-6), and the applied peRETSPL values have therefore been derived from two studies: (1) by PTB in Germany (2008), and (2) by DTU in Denmark (Gøtsche-Rasmussen et al., 2012). The mean of the values obtained in the two studies, rounded to the nearest 0.5 dB, are used by Interacoustics A/S to calibrate the broad-band CE Chirp® LS and the four frequency specific NB CE-Chirp® LSs delivered by the EAR-3A earphone.



| Toneburst<br>ECochG/ABR15/ABR30/AMLR/RATE STUDY/VEMP 0 dB 2-2-2 cycle<br>linear envelope |              |           |      | Toneburst<br>ALR/MMN dB 25-50-25 ms            |              |           |      |
|--|--------------|-----------|------|--|--------------|-----------|------|
| Hz   | Insert phone | Headphone | Bone | Hz   | Insert phone | Headphone | Bone |
| 250  | 28.0         | 38.0      | 74.5 | 250  | 17.5         | 27.0      | 67.0 |
| 500  | 23.5         | 25.5      | 69.5 | 500  | 9.5          | 13.5      | 58.0 |
| 750  | 21.0         | 23.0      | 61.0 | 750  | 6.0          | 9.0       | 48.5 |
| 1000   | 21.5         | 21.5      | 56.0 | 1000   | 5.5          | 7.5       | 42.5 |
| 1500   | 26.0         | 23.0      | 51.5 | 1500   | 9.5          | 7.5       | 36.5 |
| 2000   | 28.5         | 24.5      | 47.5 | 2000   | 11.5         | 9.0       | 31.0 |
| 3000   | 30.0         | 26.5      | 46.0 | 3000   | 13.0         | 11.5      | 30.0 |
| 4000   | 32.5         | 32.0      | 52.0 | 4000   | 15.0         | 12.0      | 35.5 |
| 6000   | 36.5         | 37.5      | 60.0 | 6000   | 16.0         | 16.0      | 40.0 |
| 8000   | 41.0         | 41.5      | 65.5 | 8000   | 15.5         | 15.5      | 40.0 |
| ISO 389-6:2007   |              |           |      | ISO 389-1:2000, ISO 389-2:1994, ISO 389-3:1994 |              |           |      |
| Click<br>ECochG/ABR15/ABR30/AMLR/RATE STUDY/VEMP 0 dB                                    |              |           |      | Click<br>ALR/MMN 0 dB                          |              |           |      |
|  | Insert phone | Headphone | Bone |  | Insert phone | Headphone | Bone |
| Click  | 35.5         | 30.0      | 51.5 | Click  | 35.5         | 30.0      | 51.5 |
| NB CE-Chirp® LS<br>ECochG/ABR15/ABR30/AMLR/RATE STUDY/VEMP 0 dB                          |              |           |      | NB CE-Chirp® LS<br>ALR/MMN 0 dB                |              |           |      |
| Hz   | Insert phone | Headphone | Bone | Hz   | Insert phone | Headphone | Bone |
| 500  | 25.5         | 27.5      | 74.0 | 500  | 25.5         | 27.5      | 74.0 |
| 1000   | 24.0         | 24.0      | 61.0 | 1000   | 24.0         | 24.0      | 61.0 |
| 2000   | 30.5         | 26.5      | 50.0 | 2000   | 30.5         | 26.5      | 50.0 |
| 4000   | 34.5         | 34.0      | 55.0 | 4000   | 34.5         | 34.0      | 55.0 |
| CE-Chirp® LS<br>ECochG/ABR15/ABR30/AMLR/RATE STUDY/VEMP 0 dB                             |              |           |      | CE-Chirp® LS<br>ALR/MMN 0 dB                   |              |           |      |
|  | Insert phone | Headphone | Bone |  | Insert phone | Headphone | Bone |
|  | 31.5         | 26.5      | 51.0 |  | 31.5         | 26.5      | 51.0 |

Only tone burst correction values change for ALR & MMN testing. For Click and CE-Chirps® LS, the same correction is applied.

Table 1



### Stimulus maximum intensity

From software 4.5, the stimulus maximal is increased for all transducers. Insert earphones and headphone can now go even louder.

To get the increased intensity maximal for the bone conductor, first, take the B81 BC, then ensure that the correct bone vibrator is chosen in the calibration setup (this allows more output for the bone vibrator).

When using a new transducer, always make sure that it is calibrated prior to use, and follow the procedure described in this service manual.

The below table is an overview of the minimum intensity performance of the various transducers when using software 4.5. Individual systems may perform even louder as this depends on the individual transducer sensitivity per frequency.

| Stimulus          |      | ABR3A       |      | DD45<br>short 2-1- |      | DD45S<br>short 2-1- |      | B81         |      |
|-------------------|------|-------------|------|--------------------|------|---------------------|------|-------------|------|
|                   |      | short 2-1-2 | Long | 2                  | Long | 2                   | Long | short 2-1-2 | Long |
| Burst             | 250  | 105         | 115  | 105                | 110  | 105                 | 110  | 50          | 55   |
| Burst             | 500  | 110         | 120  | 115                | 120  | 115                 | 120  | 70          | 80   |
| Burst             | 750  | 110         | 120  | 120                | 120  | 120                 | 120  | 70          | 85   |
| Burst             | 1000 | 110         | 120  | 120                | 120  | 120                 | 120  | 75          | 90   |
| Burst             | 1500 | 110         | 120  | 115                | 120  | 115                 | 120  | 80          | 95   |
| Burst             | 2000 | 110         | 120  | 115                | 120  | 110                 | 120  | 75          | 90   |
| Burst             | 3000 | 110         | 120  | 120                | 120  | 115                 | 120  | 65          | 85   |
| Burst             | 4000 | 105         | 120  | 115                | 120  | 110                 | 120  | 65          | 80   |
| Burst             | 6000 | 90          | 110  | 100                | 120  | 100                 | 120  | 45          | 65   |
| Burst             | 8000 | 70          | 95   | 95                 | 120  | 90                  | 115  | 35          | 60   |
| CE-Chirp          |      | 100         | 105  | 110                | 110  | 110                 | 110  | 70          | 70   |
| Click             |      | 100         | 100  | 105                | 105  | 105                 | 105  | 70          | 70   |
| Click 200Hz-10kHz |      | 95          | 95   | 105                | 105  | 105                 | 105  | 70          | 70   |
| NB CE-Chirp       | 500  | 105         | 105  | 115                | 115  | 115                 | 115  | 60          | 60   |
| NB CE-Chirp       | 1000 | 110         | 110  | 115                | 115  | 120                 | 120  | 70          | 70   |
| NB CE-Chirp       | 2000 | 105         | 105  | 115                | 115  | 110                 | 110  | 70          | 70   |
| NB CE-Chirp       | 4000 | 105         | 105  | 115                | 115  | 110                 | 110  | 60          | 60   |



## TEOAE

|                         |  |  |
|-------------------------|--|--|
| <b>Medical CE-mark:</b> | The CE-mark in combination with MD symbol indicates that Interacoustics A/S meets the requirements of the Medical Device Regulation (EU) 2017/745 Annex I<br>Approval of the quality system is made by TÜV – identification no. 0123 |  |
| <b>Standards:</b>       | <b>Test signal</b>   | IEC 60645-3:2007   |
|                         | <b>OAE</b>   | TEOAE IEC 60645-6:2022, Type 1 & 2 Otoacoustic emissions   |
| <b>Stimulus:</b>        | <b>Type:</b>   | Click Non-linear   |
|                         | <b>Bandwidth:</b>  | 500 – 5500 Hz  |
|                         | <b>Level:</b>  | 30 to 90 dB peSPL, peak to peak calibrated, AGC controlled |
|                         | <b>Level Step:</b>   | 1 dB SPL   |
|                         | <b>Transducer:</b>   | Dedicated DPOAE/TEOAE probe (Accuracy 0.5 dB)              |
| <b>Recording:</b>       | <b>Analysis time:</b>  | 5 seconds to 30 minutes                                    |
|                         | <b>Sampling frequency</b>  | 30 kHz   |
|                         | <b>A/D Resolution:</b>   | 16 bit, 3.7 Hz resolution                                  |
|                         | <b>Artifact Reject System:</b>   | 0 to +60 dB SPL or off Applicable during testing           |
|                         | <b>SNR Criteria:</b>   | Adjustable between 5 and 25 dB                             |
| <b>Display gain:</b>    | <b>General display gain:</b>   | Applicable during testing                                  |

|                                  |                     |  |
|----------------------------------|---------------------|--|
| <b>OAE Probe Specifications:</b> |                     |  |
| <b>Probe:</b>                    | <b>Application:</b> | TEOAE measurements   |
|                                  | <b>Dimensions:</b>  | (W x D x H) 12 x 26 x 11 mm (exc. Eclipse)                         |
|                                  | <b>Weight:</b>      | 3 g (exc. Cable, exc. Eclipse)<br>39 g (incl. cable, exc. Eclipse) |
| <b>Cable:</b>                    | <b>Length:</b>      | 2980 mm cable  |

### NOTICE

The OAE probe settings match the individual Eclipse, as the calibration data is stored on the instrument. Using other transducers require individual calibrations of the transducers.

**TEOAE**  
Probe stimuli are calibrated in peSPL values using an ear simulator coupler made in accordance with IEC 60318-4



## DPOAE

|                         |  |   |
|-------------------------|--|---|
| <b>Medical CE-mark:</b> | The CE-mark in combination with MD symbol indicates that Interacoustics A/S meets the requirements of the Medical Device Regulation (EU) 2017/745 Annex I<br>Approval of the quality system is made by TÜV – identification no. 0123 |   |
| <b>Standards</b>        | <b>Test Signal:</b>  | IEC 60645-1:2012 /ANSI S3.6   |
|                         | <b>OAE</b>   | DPOAE IEC 60645-6:2022, Type 2 Otoacoustic emissions  |
| <b>Stimulus:</b>        | <b>Frequency Range:</b>  | 500-10000 Hz  |
|                         | <b>Frequency Step:</b>   | 25 Hz   |
|                         | <b>Level:</b>  | 30 to 70 dB SPL   |
|                         | <b>Level Step:</b>   | 1 dB SPL  |
|                         | <b>Transducer:</b>   | Dedicated DPOAE/TEOAE probe   |
| <b>Recording:</b>       | <b>Analysis time:</b>  | minimum 2 sec to unlimited test time  |
|                         | <b>A/D Resolution:</b>   | 16 bit, 3.7 Hz resolution   |
|                         | <b>Sampling frequency</b>  | 30 kHz  |
|                         | <b>Artifact Reject System:</b>   | -30 to +30 dB SPL or off. Applicable during testing   |
|                         | <b>Stimulus Tolerance:</b>   | Adjustable between 1 and 10 dB  |
|                         | <b>SNR Criteria:</b>   | Adjustable between 3 and 25 dB  |
|                         | <b>Probe check window</b>  | 256 points frequency response of the ear canal due to a click stimulus presented with a rate of 100 Hz at 80 dB SPL |
|                         | <b>DP-Response window</b>  | 4096 points frequency response  |
| <b>Display gain:</b>    |  |   |
|                         | <b>General display gain:</b>   | Applicable during testing   |

|                                  |                     |  |
|----------------------------------|---------------------|--|
| <b>OAE Probe Specifications:</b> |                     |  |
| <b>Probe:</b>                    | <b>Application:</b> | DPOAE measurements   |
|                                  | <b>Dimensions:</b>  | (W x D x H) 12 x 26 x 11 mm (exc. Eclipse)                         |
|                                  | <b>Weight:</b>      | 3 g (exc. Cable, exc. Eclipse)<br>39 g (incl. cable, exc. Eclipse) |
| <b>Cable:</b>                    | <b>Length:</b>      | 2980 mm cable  |

### NOTICE

The OAE probe settings match the individual Eclipse, as the calibration data is stored on the instrument. Using other transducers require individual calibrations of the transducers.

The DPOAE module uses an improved method of stimuli level control, which more accurately delivers the specified intensity in the full range of ear canals, from infants to adults. The applicability of the IEC 60645-6 standard is currently limited to adult ears. Therefore, in order to better serve a market with a product that provides more accurate stimulus levels to a wide range of ear canal volumes (specifically infants), we have chosen utilization of a more comprehensive calibration procedure for DPOAEs that is outside the scope of IEC 60645-6 for some protocols. This improved method of stimulus control is enabled when the *Use Microphone compensation* checkbox is checked. To apply the IEC60645-6 calibration method, uncheck the *Use Microphone compensation* in the *Advanced* tab of the protocol setup.

**DPOAE**  
Probe stimuli L1 and L2 are calibrated individually in SPL values using an ear simulator coupler made in accordance with IEC 60318-4



## ABRIS

|  |  |   |
|--|--|---|
| <b>Medical CE-mark:</b>                              | The CE-mark in combination with MD symbol indicates that Interacoustics A/S meets the requirements of the Medical Device Regulation (EU) 2017/745 Annex I<br>Approval of the quality system is made by TÜV – identification no. 0123 |   |
| <b>Standards:</b>                                    | <b>Test signal</b>   | EC 60645-3:2007   |
|  | <b>AEP</b>   | IEC 60645-7:2009 Type 2   |
| <b>EPA Preamplifier:</b>                             | <b>Two channels standard:</b>  | EPA4 Cable Collector (4 electrodes). Standard 50 cm. Option: 5 cm or 290 cm                 |
|  | <b>One Channel (optional):</b>   | EPA3 Cable Collector (3 electrodes). 50 cm  |
|  | <b>Gain:</b>   | 80 dB   |
|  | <b>Frequency response:</b>   | 0,5 – 11.3 kHz  |
|  | <b>CMR Ratio:</b>  | Minimum 100 dB. Typical 120 dB @55 Hz   |
|  | <b>Noise (RTI)</b>   | =< 15 nV/√Hz  |
|  | <b>Radio frequency immunity:</b>   | Typically 20 dB improvement over previous available designs                                 |
|  | <b>Max input offset voltage:</b>   | 2,5 V   |
|  | <b>Input impedance:</b>  | >=10 MΩ/ =<170 pF   |
|  | <b>Power from main unit:</b>   | Insulated power supply with 1500 V isolation. The signal is digitally/capacitive insulated. |
| <b>Specifications as EPA4 Impedance measurement:</b> |  | Selectable for each electrode   |
|  | <b>Measurement frequency:</b>  | 33 Hz   |
|  | <b>Waveform:</b>   | Rectangular   |
|  | <b>Measurement current:</b>  | 19µA  |
|  | <b>Range:</b>  | 0.5 kΩ – 25 kΩ  |
| <b>Stimulus:</b>                                     | <b>Stimulus rate:</b>  | 93 Hz   |
|  | <b>Level:</b>  | 30, 35, 40 dBnHL  |
|  | <b>Click:</b>  | 100 µs  |
| <b>Recording:</b>                                    | <b>Analysis time:</b>  | 120 seconds   |
|  | <b>A/D resolution:</b>   | 16 bit  |
|  | <b>Sampling frequency</b>  | 30 kHz  |
|  | <b>Artifact rejection system:</b>  | Standard voltage based system   |
| <b>Display:</b>                                      | Stimulus level and type, Graph view  |   |
| <b>Security:</b>                                     | Password protection of test parameters possible.   |   |
| <b>Algorithmic Sensitivity:</b>                      | <b>Click:</b>  | 99.99%  |
| <b>Specificity:</b>                                  | <b>Click:</b>  | ≥ 97%   |



## ASSR


|                                 |   |   |
|---------------------------------|---|---|
| <b>Medical CE-mark:</b>         | The CE-mark in combination with MD symbol indicates that Interacoustics A/S meets the requirements of the Medical Device Regulation (EU) 2017/745 Annex I<br>Approval of the quality system is made by TÜV – identification no. 0123. |   |
| <b>Standards:</b>               | <b>Test signal:</b>   | IEC 60645-3:2007  |
|                                 | <b>AEP</b>  | IEC 60645-7:2009, Type 1.   |
| <b>EPA Preamplifier:</b>        | <b>Two channels standard:</b>   | EPA4 Cable Collector (4 electrodes). Standard 50 cm. Option: 5 cm or 290 cm                                       |
|                                 | <b>One Channel (optional):</b>  | EPA3 Cable Collector (3 electrodes). 50 cm  |
|                                 | <b>Gain:</b>  | 80 dB   |
|                                 | <b>Frequency response:</b>  | 0,5 – 11.3 kHz  |
|                                 | <b>CMR Ratio:</b>   | Minimum 100 dB. Typical 120 dB @55 Hz   |
|                                 | <b>Noise (RTI)</b>  | =< 15 nV/√Hz  |
|                                 | <b>Radio frequency immunity:</b>  | Typically 20 dB improvement over previous available designs   |
|                                 | <b>Max input offset voltage:</b>  | 2,5 V   |
|                                 | <b>Input impedance:</b>   | 10 MΩ/ 170 pF   |
| <b>Impedance measurement:</b>   | <b>Waveform:</b>  | Rectangular   |
|                                 | <b>Measurement current:</b>   | 19µA  |
|                                 | <b>Range:</b>   | 0.5 kΩ – 25 kΩ  |
| <b>Stimulus:</b>                | <b>Stimulus rate:</b>   | 40 or 90 Hz   |
|                                 | <b>Transducer:</b>  | Ear Tone ABR insert phone, calibrated on an IEC 711 coupler.<br>Headphone (optional)<br>Bone conductor (optional) |
|                                 | <b>Level:</b>   | 0 – 100 dB nHL in 5 dB steps.   |
|                                 | <b>NB CE-Chirp® Freq.:</b>  | 500, 1000, 2000, and 4000 Hz, both ears same time.  |
|                                 | <b>Bandwidth:</b>   | 1 octave ± ½ octave – 3 dB  |
|                                 | <b>Masking:</b>   | White noise 0 – 100 dB SPL  |
|                                 | <b>Analysis Time:</b>   | 6 minutes to detect a ASSR signal – can be extended up to 15 minutes  |
|                                 | <b>Recording:</b>   | <b>Sampling frequency:</b>  |
| <b>Artifact Reject System:</b>  |   | Standard voltage based system   |
| <b>Gain:</b>                    |   | 74 – 110 dB. Auto or Manual selection.  |
| <b>Channels:</b>                |   | 2, with separate detection algorithm  |
| <b>Algorithmic Sensitivity:</b> |   | 99% or 95% , false pass probability   |
| <b>Rejection levels:</b>        |   | Manual 5, 10, 20, 40, 80, 160, 320, 640 µV input  |
| <b>Anti- aliasing filter:</b>   |   | Analog 5kHz, 24 dB / octave   |
| <b>Display:</b>                 | Independent control of up to 8 simultaneous stimuli (max 4 per ear)   |   |
| <b>Display Gain:</b>            | Independent start, stop control for each of the 8 stimuli   |   |
| <b>Controlled parameters:</b>   | Stimulus level control for each of the 8 stimuli  |   |
|                                 | False pass probability 1 or 5%  |   |
|                                 | Test protocols included for children and adult  |   |
| <b>NOAH:</b>                    | NOAH 4 compatible   |   |

## NOTICE

The transducer settings match the individual Eclipse, as the calibration data is stored on the instrument. Using other transducers require individual calibrations of the transducers.



## Appendix B - exploded view




- 1. Upper cabinet
- 2. PCA mainboard
- 3. Diode
- 4. Lower cabinet
- 5. Power supply
- 6. Rear panel












## Appendix C - spare parts & accessories/disposables

Exploded view of Spare parts:
















**When ordering spare parts, please supply SN of the product to ensure the correct parts to be supplied.**

| POSITION NO. | SPARE PARTS   | PART NUMBER | DESCRIPTION                 |
|--------------|---|-------------|-----------------------------|
| 1            |    | 8001498     | Upper cabinet               |
| 2            |    | 8503874     | Rear panel                  |
| 3            |    | 8006026     | Power supply                |
| 4            |    | 8001239     | Guard plate                 |
| 5            |  | 8009908     | Power switch                |
| 6            |  | 8514201     | PCA mainboard 1.1 1086      |
| 7            |  | 8001461     | Lower cabinet               |
| 8            |  | 8001790     | Bumpers                     |
|              |  | 8104639     | Ground cable – internal     |
| 12           |   | 8003041     | M3x8 torx pan head          |
| 13           |   | 8003218     | M3x2.4 hex nut              |
| 18           |   | 8003285     | Ø3.2xØ7x0.5 flat w. chamfer |












| ACCESSORIES<br>DISPOSABLES /  | PART NUMBER | DESCRIPTION   |
|---|-------------|---|
|    | 8011175     | OAE probe   |
|    | 8101892     | IP30 insert phone ABR                                   |
|    | 8107307     | DD45 shielded headset with HB7 headband                 |
|   | 8107307     | DD45 shielded headset with HB7 headband – non-stretched |
|  | 8517075     | B81 bone conductor headset                              |
|  | 8103467     | Accessory kit EP  |
|  | 8103992     | Clothing clip kit                                       |
|  | 8513986     | eABR trigger kit Eclipse 1.1                            |











| ACCESSORIES<br>DISPOSABLES  | PART NUMBER | DESCRIPTION   |
|---|-------------|---|
|    | 8013229     | EPA3 cable collector 500 mm                                 |
|    | 8013232     | EPA4 cable collector 2.9 m                                  |
|    | 8013231     | EPA4 cable collector 50 mm                                  |
|   | 8507803     | Shielded preamplifier EPA4 cable connector 500 mm           |
|  | 8507804     | Shielded preamplifier EPA4 cable connector 500 mm and 50 mm |
|  | 8507805     | Shielded preamplifier EPA3 cable connector 500 mm           |
|  | 8532394     | Standard preamplifier EPA4 cable connector 500 mm           |
|  | 8533838     | Standard preamplifier EPA4 cable connector 500 mm and 50 mm |











| ACCESSORIES /<br>DISPOSABLES   | PART NUMBER | DESCRIPTION   |
|--|-------------|---|
|   | 8533841     | Standard preamplifier EPA3 cable connector 500 mm                     |
|   | 8011241     | Cable USB 2m  |
|   | 8011252     | LBK15 loop back   |
| <br> |             | Eartips of various sizes<br>Types: umbrella, mushroom, flanged        |
|   | 8013209     | TM electrode for EChoG starter kit                                    |
|   | 8101622     | IA mouse pad  |
|   | 8011348     | Probe cleaning tool<br>Assembled / open                               |
|   | 8011380     | Assortment box BET 25<br>TEOAE/DPOAE<br>Including probe cleaning tool |








| ACCESSORIES /<br>DISPOSABLES  | PART NUMBER | DESCRIPTION  |
|---|-------------|--|
|    | 8500140     | Nipple tube kit (for earphones), 20 pcs                        |
|    | 8012960     | Probe tip kit (for OAE probe, Ø3,3), 8 pcs                     |
|    | 8500370     | Gold tip trode, 13 mm, 20 pcs                                  |
|   | 8500440     | Gold tip trode, 10 mm  |
|  | 8507327     | Blue and red replacement silicone tube kit for IP30™ (Sanibel) |
|  | 8537295     | EARTone replacement velcro col                                 |
|  | 8500330     | Ten-20 electrode paste, conductive gel                         |
|  | 8500320     | NuPrep gel 4 oz/114 g tube, preparation gel                    |



| ACCESSORIES<br>DISPOSABLES  | PART NUMBER | DESCRIPTION                      |
|---|-------------|----------------------------------|
|    | 8500050     | Alcohol pads                     |
|    | 8500040     | Gauze swabs                      |
|    | 8500740     | Spectra 360 Electrode gel        |
|    | 8012952     | Abrasive paste, preparation gel  |
|  | 8003566     | Floss bridge implant             |
|  | 8500380     | Sonavelle gel                    |
|  | 8107483     | SP90A speaker incl accessories   |
|  | 8600546     | Adapter 3.5 to 6.5mm stereo jack |



| ACCESSORIES<br>DISPOSABLES /  | PART NUMBER | DESCRIPTION   |
|---|-------------|---|
|    | 8533283     | Speaker stand SP90A                                       |
|   | 8536952     | Sound field microphone                                    |
|  | 8534403     | Microphone stand ( <b>without</b> microphone adapter)     |
|  | 8603923     | Microphone adapter  |
|  | 8527167     | Loudspeaker cable 10m (6.3mm jack plug & 3.5mm jack plug) |



## Appendix D - update news

The following modifications have been made to the instrument and/or the manual:

| DATE    | ACTION  | REMARKS  |
|---------|---|--|
| 2025/07 | Revision of service manual  | - Appendix C updated with new part number from 8600014 to 8603923.   |
| 2024/05 | Revision of service manual  | - Section 4.4 Aided Cortical calibration added<br>- Appendix A: Technical specifications updated with Aided Cortical<br>- Appendix C: part numbers updated                               |
| 2024/02 | Revision of service manual  | - Appendix A: Technical specifications updated   |
| 2023/02 | Revision of service manual  | - Block diagram removed<br>- ABRIS added under 4.2<br>- 8500720 replaced by 8507327  |
| 2022/06 | Revision of service manual  | - Part numbers are updated<br>- Fax number is removed  |
| 2022/01 | Spare part list is updated and aligned against BOM. Exploded view of spare part list is included.   |  |
| 2020/02 | Revision of service manual for Eclipse (8507843): This is a part of the Rabbit software update project: IA OAE Suite 1.1 - new Eclipse FW<br><br>This revision concerns:<br>- Appendix A: New tech specs added.<br>- 4.2.2: Addition of research stimuli content - provided by KLIN (no relation to Rabbit)<br>- 3.4: OtoAccess® Database start-up updated to new format.<br>- smaller visual updates<br>- Appendix C: Addition of LBK loop back: 8011252<br>Appendix C: Correction of PCA mainboard part no.: 8513815 -> 8514201   | Valid from software version EP15/25 v4.6, IA OAE suite 1.1, ABRIS v1.06.1, VEMP v4.6 and ASSR v 1.2.8.<br><br>PCA 8514201: From serial SN0958585 valid with first release of this manual |
| 2019/06 | Revision Eclipse 1.1 - service manual:<br>- Addition of OAE probe in calibration chapter.<br>- Update to revised technical specifications.<br>- Parts & Accessories list:<br>- EABR kit: new part number (not yet provided by PM), and new photo to display mini-jack.<br>- Addition of EABR kit with cable (no part number created yet.<br>- Assortment box: revision: 8012838 ->8011380 + new photo.<br>- Mainboard: revision to calibrated mainboard kit, from previous 8108687 to 8513815.<br>- New screenshot: 0 dB should be added.<br><br>Addition of ASSR calibration paragraph.<br>Addition of Stimulus Maximum Intensity passage in appendix A.<br>New CVI. |  |
| 2017/12 | Release of service manual   |  |