Instructions for Use

Part 1

Instructions for Use – GB
Manual de instrucciones – ES
Bruksanvisning – SE
Instructions d'utilisation – FR
Οδηγίες Χρήσης - GR
Gebrauchsanweisung – DE
Istruzioni per l'uso – IT
Instruções de utilização – PT

VNG Stand Alone Unit VN415 / VO425

Valid from software version 7.07
8011826 Rev. 3 - 11/2013

Instructions for Use - GB

VNG Stand Alone Unit VN415 / VO425

Valid from software version 7.00

Intended Use

The two Interacoustics VNG systems, VN415 and VO425, are primarily intended for use in the evaluation, documentation and diagnosis of vestibular disorders. This is of particular interest to Ear, Nose and Throat specialists, to Neurology specialists, Audiology specialists and other health care professions concerned with balance assessment. The VN415 is the basic unit for peripheral vestibular testing and the VO425 is the extended system for central testing using oculomotor stimuli.

Precautions

 WARNING	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
 CAUTION	CAUTION , used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
 NOTICE	NOTICE is used to address practices not related to personal injury.

NOTICE

1. It is recommended that parts which are in direct contact with the patient (e.g. disposable goggle foam pads) should only be used with one patient, and should therefore be discarded after each session.
2. Please note that the CE Marking is only legal if this Instruction is translated into the national language of the user no later than at the delivery to him, if the national legislation demands a text in the national language according to MDD article 4.4.
3. Please note that if connection is made to standard equipment like printers and network, special precautions must be taken in order to maintain medical safety.
4. Although the instrument fulfils the relevant EMC requirements precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g. from mobile phones etc. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears.

If this apparatus is connected to one or more other devices with medical CE marking, to make up a system or pack, the CE marking is only valid also for the combination if the supplier has issued a declaration stating that the requirements in the Medical Device Directive article 12 are fulfilled for the combination.

Operation Instructions

Preparations prior to the VNG Test

To ensure reliable and valid test results from your VN415 / VO425 it is important that the patient is prepared for the test.

Preparation of the patient

No glasses or contact lenses should be worn during the examination. Eye makeup should also be removed, as this may disrupt the registration of the pupil. Finally, the patient should have received instructions to discontinue any form of non life sustaining medication prior to the examination.

Calibration

To make an allowance for different depths of patient eye orbits a calibration procedure can be performed.

This is initiated using the  button. The setup of this procedure can be altered in the configuration screen. Calibration points can be presented to the patient via the VisualLab software or via dots on the wall.

Placement of the mask

The mask should be placed on the head. Tighten the headband around the head.

Make sure that the pupil is visible and in the centre of the eye image. The Infra-red mirrors inside the cover have two possible positions. Make sure that they are adjusted to the most suitable position.

Adjust the cameras with the vertical and horizontal adjustment knobs. The image must be brought into focus with the focus knob – always avoid unfocused, blurred images since this will disrupt eyetracking.

Setting of the pupil identification threshold

The software will automatically detect the optimum grey scale for pupil identification, but this may also be done manually. To adjust the threshold manually simply find the optimum position for the slider controlling the setting for the grey-scale.

This is done by tweaking the  and  arrows on the sliders either side of the eye images.

Short Instruction


- **HELP**

The Operation Manual is a very important partner in getting to know your VNG unit.

This is available in hard copy and Online in the Electronic operation manual: Select Help/Help in the main menu bar, and a comprehensive list of help topics makes up an operation manual including easy searching and cross referencing functions.

- **Operating the OtoAccess™ database**

The first screen that will come up after the system has been turned on is the OtoAccess™ database screen.


The screenshot shows the OtoAccess database interface. It features a menu bar with 'File', 'Edit', 'View', 'Help', and 'Login'. Below the menu is a toolbar with icons for 'Clients' and 'Journal'. The main area is divided into several sections: 'First name' (Demo), 'Last name' (Demo), 'Initials' (DMO), 'Person ID' (0101013333), 'Birth date' (31-05-1970), 'Gender' (Male), 'Phone' (+45 6371 3555), 'Mobile' (empty), 'Fax' (+45 6371 3522), 'Email' (empty), 'Picture' (a small icon), 'Address' (Drøjerveangot 8), 'Address2' (empty), 'Zipcode' (DK-5610), 'City' (Assens), 'Country' (Denmark), 'State' (empty), 'County' (FYN), and a checkbox for 'With schedule today'. On the right side, there is a 'Select Instrument' section with a dropdown menu showing 'VO425', and a 'Categories' section with two dropdown menus for 'Category 1' and 'Category 2'.

To add a new client to the database click the “New Client” button  and after entering the required details such as name, person ID etc, save by .

To search for any client, click the “Search Client” button , write client details (e.g. Person ID) and use Enter.

To add a new examiner, open the File tab and click “Setup”  **Setup...** in the pull-down menu. Open the “Examiners” panel and enter the relevant details for the new examiner. Left click on the “Save” button to store the examiner.

• New Test

To enter test mode from the name database, double click the instrument symbol in the “Select Instrument” box .

Open a new session from the Session Template by clicking , or simply click one of the test icons to prepare for immediate testing.

Start the test with the  button


Then use the return key or the IR remote control to progress through the steps of each test. The IR remote control acts as a functional return key and can be used to start and stop testing, and for switching on the fixation light.

• Activating the fixation light

To activate the fixation light inside the goggle the  can be used. Automatic activation can be set in the Configuration dialogues.


• Manual Edit Mode

To enter manual editing mode, click the “Toggle edit mode” button . In this mode graphs are displayed in single mode, and in a larger scale than in the normal mode.

• Video recording

With the test in progress, click the ‘Start Record’ button located on the left of the eye images. To change the default setting of a test into automatic video recording, this option must be activated in the test configuration dialogue under ‘Configure Test’.

When ‘Record’ button has been activated it will flash until the test is stopped, or until you press the ‘Stop recording’ button. Recording will automatically stop when the test is stopped or finished.

To review a recorded test, activate the  button in the toolbar. The screen will now change into playback mode and the cameras will be disconnected.

- **Adding notes**

Notes about the Session or the individual test can be added by clicking on the  button.

- **Creating a Letter**

Select the 'Letter Wizard. . .' from the 'Session' menu. Select the appropriate information from the presented fields. New selectable information can be added with the  button. This can then be saved with  button or deleted with  button.

- **Printing**

Select the "Print Session" button  to receive printouts from the whole session. Select the "Print Test" button  to receive printouts for a selected test.

- **Save and Exit**

Select  to save data and  to exit session. Alternately select  if you want to exit without saving data.

Various additional tools:


Lock Chair


Click the 'Lock Chair' button  to lock the Interacoustics rotary chair in any horizontal position. This function is relevant for clinics using the Interacoustics rotary chair as an examination table for caloric irrigations and positional testing. To release the chair, click the 'Lock Chair' button again.

Tree View

Click , to access the tree view. This will display the Session and test order information on the left side of the screen.

Session Template Designer

To view and make changes to a Session Template choose 'Session Template Designer' under the Session menu. In this dialogue the properties of a Session can be changed. A new Session Template can be configured by pressing the  button. Changes to the selected template can be saved with the  button. The number of tests performed and their order can be changed in

the tests to be performed list . Tests can be inserted by using the  button. Tests can be deleted using the  button. The order of tests performed can be changed by using the  and  buttons. The settings of a selected test can be changed in by using the  button.

Automatic threshold detection

The grey-scale threshold by which the pupil is identified is set automatically when the goggle is positioned on the client. This can be deactivated and then re-activated with the  and  buttons.

Shifting between Eye Trackers

To facilitate optimal eyetracking on as many patients as possible, two different eyetrackers are available for the user:

- **Elliptical** – searches for edges; those that combine into an elliptical shape are defined as the pupil. Will normally provide optimal tracking on patients wearing mascara
- **Circular** – searches for dark areas and locks onto the roundest object, identifying this as the pupil.

To shift between the two eye trackers on the fly, click the  symbol in the lower left corner of the eye image area and select the required eye tracker – Elliptical or Circular.



Image Brightness Adjustment

Move the sliders to adjust the image brightness in enlarged eye view mode (opened by right-clicking the eye image).

The eye tracker will have optimal working conditions on images with significant contrasts between the pupil and the iris.

Change test configuration


When a test has been selected the configuration of the test can be changed using the  button. Alternately, the eye position graphs can be double clicked on. In the dialogue box that is displayed, test/graphs/analysis/stimulation/calibration parameters can be changed.

Change nystagmus identification

The thresholds by which nystagmus analysis is determined and the sensitivity of slow phase detection can be set in the configuration dialogue under analysis. If changes in these values are to be saved the  button can be used. If changes are not saved and to use the last saved changes press the  button.

Record different traces

To record different eye traces check the appropriate radio button in the test panel

section of the Configuration screen.

Changes to oculomotor stimuli (only with VO425)

Changes can be made to the Oculomotor stimuli characteristics in the stimulus leaf of the configure screen.

Testing Oculomotor Stimuli (only with VO425)

Stimuli can be previewed without recording by using the  button within the Configuration/Stimuli screen.

Changes to geometry of test environment (only with VO425)

The geometry of the test environment can be changed to accommodate different test situations. To change these settings new horizontal, vertical and distance from eye to display measurements must

be entered in the  dialogue. An outline from which to measure the vertical and


horizontal projection sizes can be given with the  button.

Connection Layout

Connect the Hardware

Connect the VNG equipment to the video projector to the laptop or to the tower PC before powering up the components.

Connect Hardware to Laptop using FireWire™ PC ExpressCard/34

1. Plug the laptop power supply cable **(1)** into the power supply and into the port at the rear of the laptop.
2. Plug the power supply cord **(2)** into a standard power outlet.
3. Plug the FireWire™ PC ExpressCard/34 **(3)** into designated slot in the side of the laptop. Do not plug it into the power supply at this time.
4. Plug the video camera cables **(4)** into the FireWire™ ports (red cable to red port; blue cable to blue port) on the FireWire™ PC ExpressCard/34.
5. Connect the power supply cable **(5)** to the FireWire™ PC ExpressCard and to a standard power outlet.
6. Plug the projector power cable **(6)** into the VGA ports on the projector and the laptop.
7. Plug the projector power cable **(7)** into a standard power outlet and into the projector power port. Do not turn on the power at this time.
8. Plug the USB connected foot pedal **(8)** into the USB port on the laptop.
9. Plug the printer cable **(9)** into a USB port on the laptop. Plug the printer power cord **(10)** into a standard power outlet.
10. Turn on the laptop, the video projector, and the printer before starting the VNG software.

NOTICE


Always power off the FireWire™ PC ExpressCard/34 before removing it from the laptop.

To connect the VNG hardware to a desktop PC, the exact same procedure as for the laptop should be followed.

In case the desktop PC does not feature an analog HD15 VGA port, or if only one DVI-I port is available, a splitter cable with the following specifications should be inserted between PC and monitor/projector:

DVI-I Male to 1 Ft. HD15 Female (VGA) / 6 Ft. DVI-D (Single Link) Male Analog and Digital Y Splitter

Connect Hardware to Tower PC

1. Plug the computer power supply cable **(1)** into the port at the rear of the computer and into the power supply.
 2. Connect the computer power supply cord **(2)** into a standard power outlet.
 3. Plug the video camera cables **(3)** into the FireWire™ ports (red cable to red port; blue cable to blue port) on the computer.
 4. Plug the main splitter cable plug **(4)** into the monitor port on the computer.
 5. Plug the monitor cable **(7)** into the splitter cable port (#1).
 6. Plug the monitor power cable **(8)** into a standard power outlet.
 7. Plug the projector cable **(5)** into the VGA ports of the projector and the splitter cable port (#2).
 8. Plug the projector power cable **(6)** into a standard power outlet and into the projector power port. DO NOT turn on the power at this time.
 9. Plug the USB connected remote control receiver **(9)** or USB connected foot pedal into a USB port on the computer.
 10. Plug the printer cable **(10)** into a USB port on the computer.
 11. Plug the printer power cord **(11)** into a standard power outlet.
 12. Turn on the computer, the video projector, and the printer before starting the VNG software.
- Note: The monitor cable must be plugged into port #1 on the splitter cable and the projector cable must be plugged into port #2 on the splitter cable.*

Explanation to symbols which can be found on the instrument:

	On (Power: connection to the mains)
	Off (Power: disconnection from the mains)
	Alternating current
	Fuse
	Ground
	Dangerous voltage
	See explanation in manual
	Type BF equipment

Technical Specifications

The technical specifications provided here cover the general aspects of the instrument, while more specific details can be found in the Operation Manual in English.

Standards

IEC 60601-1	(General Safety) Class II, Type BF
IEC 60601-1-1	(Safety of Systems)
IEC 60601-1-2	(EMC)

The CE mark indicates that Interacoustics AS meets the requirements of Annex II of the Medical Device Directive 93/42/EEC. Approval of the quality system is made by TÜV - Identification no. 0123.

Power Supply:

Input Volts:	100-240 V AC 50/60Hz 0.5A universal input switch mode.
Safety:	VDE750, EN60601-1, IEC601-1, UL544, CSA22.2

Software Features:

External patient database
Real time eye image analysis
Implemented norm data
Windows® XP operating system
Windows® 7 (32 & 64 bit)
Windows® Vista Business Premium & Ultimate

Parts included

- Interacoustics Combi Goggle w. FireWire® Camera Module(s)
- USB connected Foot Switch
- IEEE1394 PCExpressCard/34 with power supply
- Operation Manual
- VNG Quick Guide
- CE Manual
- VNG Software Installation CD
- Disposable Goggle Foam Pads, Box of 24 units

Concerning Repair

1. Interacoustics is only considered to be responsible for the validity of the CE marking, effects on safety, reliability and performance of the equipment if:

assembly operations, extensions, readjustments, modifications or repairs are carried out by authorized persons,

a 1 year service interval is maintained

the electrical installation of the relevant room complies with the appropriate requirements, and

the equipment is used by authorized personnel in accordance with the documentation supplied by Interacoustics.

2. Should a problem with the device arise, please contact your local distributor. The RETURN REPORT must also be sent to Interacoustics, Drejervaenget 8, DK-5610 Assens in connection with reporting of problems, incidents, etc.

General Information

The Rotary Chair is constructed to a maximum patient weight of 135 kilos

This instruction manual is a part of the instrument and should always be kept at hand.

It documents the state of development of the unit, valid at the time of printing, and contains all the information required by the IEC 60601-1.

Exact adherence to the instruction manual is a prerequisite for the perfect and intended functioning of the Interacoustics instrument.

Always be sure to disconnect the power supply before cleaning or disinfecting mains operated instruments.

Interacoustics A/S only guarantees for the safety, reliability and function of the instrument, if: installation, extensions, modifications and repair are exclusively carried out by personnel authorized by Interacoustics A/S;

the room in which the equipment is installed complies with the country-specific regulation;

the installation site complies with the ambient conditions of the instrument;

the unit is used in accordance with the instruction manual after initiation in the operation;

The device is not approved for operation in potentially explosive rooms or flammable atmospheres.

If fluid gets into the instrument, the instrument will have to be switched off at once. The fluid must immediately be removed by Interacoustics service technicians and the instrument tested for technical safety. Defective connectors and leads must immediately be replaced or repaired by a specialist or by personnel authorized by Interacoustics A/S.

Precautions


Following the ANSI recommendations (American National Standards Institute) for safety notes, specific passages of this instruction manual are clearly marked as safety notes.

	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
	CAUTION , used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
	NOTICE is used to address practices not related to personal injury

Operation Instructions

Examination Room Layout

The Rotary Chair may be installed using the following guidelines for room size and hardware installation:

Hardware installation

The rotary chair must be secured to the floor with 4 screws through the bottom rim of the base. Please use 12 mm dowels and 8 x 120 mm screws.

The Patient is not allowed to sit on the calves-area of the chair when this has been reclined into its horizontal position

Never connect the instrument to the power outlet before the side panels and the top cover have been correctly refitted and secured.

Consult with the product specification plate on the unit for compliance between the unit's own data and those of the local power supply system (mains voltage and mains frequency) before connecting the unit.

Connect only if all data agrees!

Inspect the mains connection cables for visible damage prior to establishing the connection. Damaged cables or plugs must immediately be replaced by Interacoustics service technicians, an authorized person or specialist.


Connection Panel

The backrest features a connection panel with purpose specific connections:

FireWire® interface for the camera(s) mounted in the goggles

USB 1.1 interface

For future use


Input/output connection panel , back rest

Rotary Chair Instruction for Use - English

Date: 2011-09-27

Page 15

Signal input/output connection panel, chair base

Software Installation

The VNG installation CD received with the VNG system holds the relevant software for the Kvaser CanBus-to-USB adapter.

The Rotary Chair can be connected to any laptop or desktop PC with the following specifications:

- Intel Core 2 Duo processor 1.8 GHz or better
- Minimum 1 GB of RAM
- Laptop PC: One 34 mm PCExpressCard slot available
- Desktop PC: Texas Instruments chipset PCI FireWire® board with dual 6-pM ports
- Three USB ports available
- One VGA port available
- Monitor resolution 1024 x 768 or better
- WindowsXP™, Windows Vista® or Windows 7® Ultimate Operating system

The Sinusoidal Pendular Test, or Sinusoidal Harmonic Acceleration (SHA)

In the Sinusoidal Pendular Rotary Chair Test the patient is positioned on a motor driven rotary chair that rotates in a sinusoidal pattern from left to right.


Figure 1: SHA test windows

Setting up the SHA test

To set up the rotary chair stimulus for the SHA test, open the 'Stimulus' tab in the configuration dialog box.

Configure Sinusoidal Pendular Test

Test | Graphs | Analysis | **Stimulus** | Calibration | Normdata

1. ☐ Manual chair control
Threshold for start: 5 [°/s]

2. ☒ Automatic chair control
Velocity [°/s]: 50 Frequency [Hz]: 0.32 Cycles: 6 Acceleration: 101 °/s²

3. ☐ Start direction left ☒ Start direction right Test time: 0:21

Velocity [°/s]	Freq [Hz]	Cycles	Test time [s]	Start direction
50	0.010	2	5:00	Right
50	0.020	2	2:30	Right
50	0.040	2	1:15	Right
50	0.080	4	1:02	Right
50	0.160	5	0:37	Right
50	0.320	6	0:21	Right

4. [New] [Update] [Remove]

8. ☒ Move chair to zero before test starts
9. ☒ Move chair to zero when test is complete
10. Zero pos: 0 [°]

[Save As Default] [Load Defaults]

[OK] [Cancel] [Help]

1. Select 'Manual chair control' if you want to control the chair manually.
2. Select 'Automatic chair control' if you have a Nydiag 200 chair connected via USB and you want VNG software to control the chair automatically.
3. Select the starting direction of the chair rotation.
4. Selected frequencies of sinusoidal harmonic acceleration testing will display in this field. The stimulus parameters will be displayed, including velocity (°/s), frequency (Hz), cycles, test time (s), and start direction.

Rotary Chair Instruction for Use - English


Date: 2011-09-27

Page 16

- Click on 'New' to add another stimulus to the test. The new stimulus will get added directly below the highlighted frequency. To make changes to the parameters of the new stimulus, highlight it in the frequencies field and enter desired velocity, frequency, and cycles in the appropriate fields.

Velocity [°/s]	Freq [Hz]	Cycles	Test time [s]	Start direction
50	0.010	2	5:00	Right
50	0.020	2	2:30	Right
50	0.020	2	2:30	Right
50	0.040	2	1:15	Right
50	0.080	4	1:02	Right
50	0.160	5	0:37	Right
50	0.320	6	0:21	Right

The minimum frequency that can be selected is .01 Hz and the maximum frequency is 2 Hz. When the desired frequency is selected, the software calculates the resulting acceleration. When the acceleration exceeds $100^\circ/\text{sec}^2$, the frequency will automatically be changed to the maximum possible.

Alternatively, by reducing the selected velocity a higher frequency can be selected. The maximum number of stimuli that can be added to one test is 10.


- Click on 'Update' for any changes made to the stimulus to become active.
- Click on 'Remove' to delete the highlighted stimulus from the list of stimuli.
- If this box is checked the chair will move to its zero position before the test is started.
- If this box is checked the chair will move to its zero position after the test is complete.
- In this field enter the number (in degrees) at which the desired zero position is for the chair.

Starting the SHA test

- To start the test click  in the test toolbar. The following dialog box will appear.

Test setup	SHA	VOR
.01 Hz - 50 °/s - 2 Cycles	☒	☐
.02 Hz - 50 °/s - 2 Cycles	☐	☐
.04 Hz - 50 °/s - 2 Cycles	☐	☐
.08 Hz - 50 °/s - 4 Cycles	☐	☐
.16 Hz - 50 °/s - 5 Cycles	☐	☐
.32 Hz - 50 °/s - 6 Cycles	☐	☐

- Select the desired SHA or VOR Suppression frequency from the list of stimuli. Once a stimuli is selected, the following dialog will display for a brief moment:

Rotary Chair Instruction for Use - English

Date: 2011-09-27

Page 17

3. If you have selected to move the chair to zero position before start of the test, the chair will automatically move to zero.
4. Once the chair is ready the above dialog will close and the test will begin.
5. Once the desired number of cycles have been completed the test will automatically stop.
6. If you have selected to move the chair to zero position after testing is complete, the chair will now move to zero.

Colour codes

- **Red lines and points** represent data of the **right eye**,
- **Blue lines and points** represent data of the **left eye**.
- **Green lines** represent data of the stimulus, in this case the **rotary chair movement**.

The Step Rotation Test

During the Step Rotation Test the patient is positioned on a motor driven rotary chair that, after an initial acceleration, rotates at a constant speed until the patient adapts to the rotation. The chair is then suddenly stopped. After a pause, the stimulus pattern is repeated, but in the opposite direction.


Figure 2: Step Rotation Test window

Setting up the Step Rotation test

To set up the Step Rotation test, open the 'Stimulus' tab in the configuration dialog box.

Rotary Chair Instruction for Use - English

Date: 2011-09-27

Page 18

1. Select 'Manual chair control' if you want to control the chair manually.
2. Select 'Automatic chair control' if you have a Nydiag 200 rotary chair connected via USB and you want VNG software to control the chair automatically.
3. When automatic chair control is selected, the Step Rotation test parameters for the highlighted stimulus will display in these fields. The minimum velocity allowed is 1°/sec and maximum is 200°/sec. The minimum acceleration/deceleration allowed is 0.1°/sec² and the maximum is 100°/sec². Stable time is the period of time in which the chair will rotate for the per-rotary test segments. Stop time is the period of time in which the chair is not rotating (post-rotary test segment). When the desired stop time has expired, the chair will automatically begin rotating the opposite direction for the next per-rotary test segment. The minimum stop time is 0 seconds and the maximum is 600 seconds.
4. Selected Step Rotation test stimuli will display here. Up to two stimuli can be added to a single Step Rotation test.
5. Click on 'New' to add another stimulus to the test. The new stimulus will get added directly below the highlighted frequency. To make changes to the parameters of the new stimulus, highlight it in the numerical graph field and enter desired velocity, acceleration, stable time, deceleration, and stop time in the appropriate fields above.

Velocity [°/s]	Accel. [°/s ²]	Stable time [s]	Deaccel. [°/s ²]	Stop time [s]
180	100	60	100	60

Test time: 4:07

Vel. [°/s]	Acc. [°/s ²]	Stable [s]	Deacc. [°/s ²]	Stop [s]	Test time [s]
50	100	60	100	60	4:02
180	100	60	100	60	4:07

New
Update
Remove

6. Click on 'Update' for any changes made to the new stimulus to become active.
7. Click on 'Remove' to delete the highlighted stimulus from the list of stimuli.
8. If this box is checked, the chair will move to its zero position before the test is started. If both right and left step rotation test segments have not yet been performed, the chair will only move to zero when the first per-rotary segment is started. The following segments will start from the chair's current position.
9. If this box is checked the chair will move to its zero position after the test is complete. If both right and left step rotation test segments have not yet been performed, the chair will only move to zero when the last test segment has been completed.
10. In this field, enter the number (in degrees) at which the desired zero position is for the chair.

Starting a Step Rotation test

1. To start a test click  in the test toolbar. The following dialog box will appear:

Select a Test

Velocity 50 °/s

☒ Rotation to the right
 ☐ Rotation to the left

OK

Velocity 180 °/s

☐ Rotation to the right
 ☐ Rotation to the left

Cancel

Save & Close


Select the desired Step Rotation velocity and rotation direction from the list.

2. Once a stimuli is selected, the following dialog will display for a brief moment:

Rotary Chair Instruction for Use - English

Date: 2011-09-27

Page 19

3. If you have selected to move the chair to zero position before start of the test, the chair will automatically move to zero.
4. Once the chair is ready the above dialog will close and test will begin.
5. Once the selected direction has been completed the test will automatically stop again. If the per-rotary segment in the opposite direction has not yet been completed the software will automatically start the per-rotary segment in the opposite direction immediately following the defined stop time for the post-rotary segment.
6. If you have selected to move the chair to zero position after testing is complete, the chair will now move to zero position.

Error Messages

The following error messages may occur during operation of the Nydiag 200 rotary chair:


Figure 0-1: Dialog appearing when no chair is connected

VNG could not detect any rotary chair. Make sure that a chair has been connected via USB and that the chair is powered on.


Figure 0-2: Patient stop button has been activated

Patient has pressed the stop button that is located under the right chair handle.


Figure 0-3: Emergency stop button has been activated

Emergency stop button has been pressed, or power to the chair has been lost.

Hygiene

Always be sure to pull the plug before cleaning or disinfecting mains operated instruments.

The instrument must be cleaned with a non-fuzzy, slightly moistened piece of cloth (please do not soak it!) Only use a mild alcoholic fluid as a disinfectant.

When cleaning the instrument, any moisture, like e.g. condensation water, must necessarily be prevented from getting inside!

Chemicals required for operation or care of the unit must always be stored, prepared and made available in specially labelled containers to prevent any mistakes.

Maintenance

Before switching on the instrument make sure that the mains cable, mains plug, socket and input at the instrument are in a perfect state. If the drive of the rotary chair becomes noisy the driving belts should be renewed. We recommend renewing the driving belts after five years.

The operator is responsible for the perfect state of the instrument. It is recommended by Interacoustics to have safety checks carried out at one-year intervals.

Immediate maintenance is required if:

**the instrument was subjected to extreme mechanical stress (impact, defective cable due to inadmissible traction);
fluid gets into the instrument;
cables and/or connectors are defective;
rubber-joints have developed cracks.**

For further details please contact your Interacoustics representative.

SYMBOLS ON INTERACOUSTICS COMPONENTS

Ground conductor connection

Unit of type B

Technical Specifications

Maximum speed	200 deg/s
Maximum acceleration	100 deg/s ²
Maximum patient weight	135 kilos
Reclining backrest	manually operated from 0 deg (horizontal) to 90 deg (upright)
Emergency stop	emergency stop button disconnects the motor power
Patient alarm button	sends an alarm signal to the computer and turns off power to the motor
Weight	195 kg
Shipping weight	230 kg
Dimensions	90 cm x 70 cm x 160 cm
Shipping dimensions	100 cm x 120 cm x 190 cm
Power supply	110-120VAC 220-240VAC
	50-60Hz 50-60Hz
	900VAC max 900VAC max

Working with ENG: EO425

Installation

During installation you need to select “ENG recording system” if using ENG as the sole recording system, and select “Firewire®/ENG recording system” if you will be using both VNG and ENG recording systems:


Figure 1.1: Set up recording system

Click OK and then select Monocular or Binocular depending on your setup. Click OK and finish the installation.

Figure 1.2 Select electrode montage

Connecting the ENG preamplifier to the computer

Connect the preamplifier via a USB cable and remember to turn the box on! The current version of the ENG preamplifier does not tell the software if it is turned on or not. It is possible to start a test with the preamplifier turned off, but you will not get any correct readings.

If the preamplifier is turned off or not connected when starting VNG you will receive this error message:


Figure 1.3 Error message when hardware is not turned on

License

Start VNG-program normally. When you start it for the first time you need to enter a license:


Figure 1.4 Error message when license not installed

Go to Tools->Setup license and you will see this dialog:


Figure 1.5 Popup box to enter serial number and license number

Enter the serial number located on the bottom of the preamplifier. If the serial number is e.g. P5590-3 you must enter 55903 in the serial number field. Do not enter characters or "-". Afterwards you enter the license. The "Store license" button will only be enabled if the license is valid.


Figure 1.6 Enter license number

Click Store license and then Exit.

Preparing patient

First you need to apply the electrodes to the patient in one of the following patterns:

2-channel left vertical


Figure 1.7 Left monocular setup

Electrode	Place
LT	Left temple
LN	Right temple
D	Below left eye
U	Above left eye
C	Forehead
RN	Shorted / Not used
RT	Shorted / Not used

2-channel right vertical


Figure 1.8: Right monocular setup

Electrode	Place
LT	Left temple
LN	Right temple
D	Below right eye
U	Above right eye
C	Forehead
RN	Shorted / Not used
RT	Shorted / Not used

3-channel left vertical


Figure 1.9: Left Binocular setup

Electrode	Place
LT	Left temple
LN	Left side of nose
D	Below left eye
U	Above left eye
C	Forehead
RN	Right side of nose
RT	Right temple

After you have applied the electrodes you need to let the electrodes “adjust to patient” for some minutes. It is different from patient to patient how much time it takes to adjust. You may experience signal drift in the first 5-10 minutes after the electrodes have been applied.

Configure electrode montage

Included in the “Configure Test” tabs for the various tests, you can select which electrode montage you will use for the testing.


Figure 1.10 Electrode montage

NOTICE

It is important to note that when switching from a monocular (2-channel) to a binocular (3-channel) electrode montage within the individual test configurations, it does not apply any electrode montage changes to the system configuration. If you wish to make a change to the default electrode montage you must make the change in the system configuration. This can be done by going to 'Tools' and then 'Switch between ENG/VNG'. A dialog box will then appear where you can select mode as either monocular or binocular.


Figure 1.11 Switch between ENG/VNG


Figure1.12 Select source

If you switch between montages, i.e. if you switch from a 2 channel montage to a 3 channel montage, you will be prompted to recalibrate, as seen in the message below:


Figure 1.13 repeat calibration

Impedance test

Once the electrodes have adjusted to the patient, you need to do an impedance test. To do this you open a test or a session, and click the calibration icon  :


Figure 1.14 Calibration and impedance check

Click on the Impedance button to start impedance test:


Figure 1.15 Impedance Check box

The white squares indicate the measured impedance. Click "Check All" to measure impedance on all electrodes. If all electrodes have impedances within tolerance, then you will have all green squares:


Figure 1.16 Impedance check complete

Green and yellow indicate acceptable. Red and black indicate unacceptable impedances. Click close to return to calibration dialog.


Figure 1.17 Electrode gain readings

The “Electrode readings” graph shows the raw readings from the electrodes. Tracings are in green, white and red as indicated on the electrode connector. When running in 2-channel mode you will only see LT/LN and D/U. If you see a lot of noise in the readings you need to check to electrodes again.

Before starting a calibration you can adjust the gain if the responses from the electrodes are too low or too high. Ask the patient to look left/right up/down to adjust the gain correctly. The tracings must not “hit” the max or min markings on the graph. If they do, you need to select a lower gain. Typically you select the medium gain setting. Filter frequency is default to 24Hz and coupling is set to DC. You do not need to change these values. Once the gain has been set you click “Start” to perform a calibration. The calibration procedure requires a minimum of 1 unit on the graph (horizontal lines) between two calibration points. So if you are having trouble to get the calibration to accept a point either your gain is too low or you have too much signal noise. You can manually click “Accept point” if the calibration cannot accept the point automatically.

It is important that you connect the electrodes correctly, do not reverse LT/LN, RN/RT, or D/U. The calibration will not work if you do so.

Calibration

After the impedance test is complete you need to calibrate. The procedure is the similar to the VNG calibration, and the patient is instructed to look at each of the 5 dots. A green check mark will appear when the calibration angle is verified.


Figure 1.18 calibration interface

Performing a test

Once the impedance and calibration have been completed there is not much difference between conducting tests in ENG or VNG. For ENG the color for left eye is blue and right is red. If you do experience signal drift while testing you can click the left foot pedal to re-center tracings. Also you can click the “Center ENG tracings” icon in the upper right corner (except for Saccades and Pursuit, which must be centered after the test is complete):


Figure 1.19 Center ENG tracings

In Smooth Pursuit and Saccade test you cannot re-center tracings manually if you experience drift. You should only make adjustments after the data collection is completed.

Smooth Pursuit drift

Below you can see an example of a smooth pursuit test with significant signal drift:


Figure 1.20 Drift in the smooth pursuit test

To correct this you need to double-click in the test window to go to Configure Smooth Pursuit Test->graph settings and enable “Compensate for baseline drift”:


Figure 1.21 Configure Smooth Pursuit Graph tab

This will adjust the tracings and re-center them:


Figure 1.22 Smooth Pursuit after compensation for drift

Saccade Test Drift

Drifting in saccade is only a “visual” problem, as the analysis is not affected by it. Below you can see a saccade test with significant signal drift. (The signal started out a zero position).


Figure 1.23 Drift in the saccade test

To re-center the tracings double-click in the saccade test window to go to Configure Saccade Test->graph settings and select “Compensate for baseline drift”:


Figure 1.24 Configure Saccade Test Graph Tab


Figure 1.25 Compensating for drift in saccade test

This will adjust the tracings and re-center them:
The longer the test is run the more likely it is that you will experience signal drift.

Accommodating for Noisy Recordings

In certain situations, it is possible that there may be noise or artifact in your ENG recordings that the software may misinterpret and mark as beats of nystagmus.


Figure 1.26 Noise in ENG recording

To accommodate for this noise so that fewer nystagmus beats are incorrectly marked, you can go into configuration for the specific test by double-clicking in the test window. Go to the Analysis tab and increase the minimum amplitude of fast phase to 2 or 3 degrees.


Figure 1.27 Configure Test Analysis Tab

With this change, the minimum amplitude of the fast phases of nystagmus will now have to be higher before the software will detect and mark beats of nystagmus. It is recommended that this minimum be left at the default setting of 1° when performing VNG recordings. You can also adjust the minimum time of the slow phase that is required for a beat to be detected. 100ms - 200ms is appropriate if you experience noise tracings.


Figure 1.28 Same recording after increasing the minimum amplitude of fast phase

Switching between ENG and FireWire® VNG

To switch between ENG and Firewire® VNG, click on Tools from the top menu bar and select 'Switch between ENG/VNG'.



Figure 9129 Switch between ENG/VNG

This may also be accomplished by exiting VNG and going to System Configuration under the Masks tab. Select Firewire® in the Data Source-section and click OK.


Figure 1.30 System Configuration Masks tab

To switch back to ENG you need to select ENG instead and click OK.

Troubleshooting

If the ENG preamplifier becomes unplugged while the software is still running you will receive the following error message:


Figure 1.31 VNG executable error message

Using the Digital Light Bar

First set up the distance from light bar to patient. Default is 1 meter or 39 inches.

The stimuli setup dialog can e.g. be accessed from the calibration dialog:


Figure 2.1 Set up DLB distance to patient

There are two different angles to consider. The first one is for saccade, smooth pursuit and OKN. These tests use the center diodes which can show stimuli at any position within the default ± 20 degrees.

The second one is used in gaze test only, as they are fixed at 30 or 25 degrees. You can only achieve 25 or 30 degrees at 1m (39") when using the long version of the light bar.

Calibration

There is only one calibration pattern available:


Figure 2.2 Calibration angles

The order of the calibration points are: 1: center, 2: left, 3: right, 4: center, 5: up, 6: down.

Click 'Start' and the calibration will calibrate point 1, 2 and 3. Then you are asked to flip the light bar to vertical position:


Figure 2.3 Vertical calibration message

Click 'OK' and point 4, 5 and 6 is completed.

Gaze test

The light bar has a default range of ± 20 degrees when using the outermost diodes. However these can only be used in gaze test. The user has to remember to flip the light bar when running the vertical stimuli.


Figure 2.4 positioning of light bar

The stimuli setup dialog now allows the user to also setup light bar. The functionality for projector is disabled when using the light bar.


Figure 2.5 Configure gaze test for use with lightbar

ENG/DLB Instruction for Use - English

Date: 2011-09-27

Page 38

Input	Description	Value range
Diodes to turn on	This allows the user to setup how many diodes to show at the time when running the test. The more diodes the easier the stimuli is to see. Default is 1.	1-4

Smooth pursuit test

This test uses the inner range of diodes which gives a default range of ± 20 degrees.

Stimuli setup dialog is identical to gaze test.

User has to remember to flip the light bar when running the vertical stimuli.

Saccade test

This test uses the inner range of diodes which gives a default range of ± 20 degrees.

Stimuli setup dialog is identical to gaze test. It is not possible to run a "mixed pattern"-test using the light bar, as it can only project stimuli in one direction at a time. User has to remember to flip the light bar when running the vertical stimuli.

Optokinetic test

This test uses the inner range of diodes which gives a default range of ± 20 degrees.

Stimuli setup dialog allow the user to modify this:


Figure 2.6 Configure test for OPK and DLB

Input	Description	Value range
Diodes to turn on	This allows the user to setup how many diodes to show at a time. E.g. when selecting 3, then the light bar turns on 3 diodes equally distanced across the light bar. Default is 3.	3-6

Light bar service tool

Go to Start menu->All programs->Interacoustics->VNG and select "DLB Service Tool" to start the service program. Make sure VNG is closed before doing this.


Figure 2.7 DLB service Tool

This will start the following program which allows you to test if all diodes work as expected.


Figure 2.8 DLB service Tool diode check

Press "Connect" to start test testing:

Figure 2.9 DLB diode functionality check

ENG/DLB Instruction for Use - English

Date: 2011-09-27

Page 40

Text	Functionality	Shortcut
Connect	Connects to the light bar.	
Close	Closes the program.	
<	Jumps to the left most diode.	"End" key
<<	Starts a sweep to the left. Press the button again to increase speed.	"Page down" key
<	Shows the next available diode to the left. Stops any sweep running.	"Left" key
>	Shows the next available diode to the right. Stops any sweep running.	"Right" key.
>>	Start a sweep to the right. Press the button again to increase speed.	"Page up" key
>	Jumps to the right most diode.	"Home" key

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
The screenshot shows the 'Demo Demo' window of the OtoAccess™ software. It features a menu bar (File, Edit, View, Help, Login) and a toolbar. The main area is divided into two panes: 'Clients' and 'Journal'. The 'Clients' pane contains a form for patient data, including fields for First name, Last name, Initials, Person ID, Birth date, Gender, Phone, Mobile, Fax, Email, Address, Address2, Zipcode, City, Country, State, and County. There is also a 'Picture' field with a thumbnail image. The 'Journal' pane shows a 'Select Instrument' dropdown menu with 'VO425' selected, and a 'Categories' section with 'Category 1' and 'Category 2' dropdowns. A checkbox labeled 'With schedule today' is at the bottom.

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The screenshot shows the 'Configure Sinusoidal Pendular Test' dialog box, specifically the 'Stimulus' tab. The dialog has several sections for configuring the test parameters.

- 1** ☐ Manual chair control
 - Threshold for start: [°/s]
- 2** ☒ Automatic chair control
 - Velocity [°/s]: Frequency [Hz]: Cycles: Acceleration: 101 [°/s²]
- 3** ☐ Start direction left ☒ Start direction right Test time: 0:21
- 4** Table with columns: Velocity [°/s], Freq [Hz], Cycles, Test time [s], Start direction.

Velocity [°/s]	Freq [Hz]	Cycles	Test time [s]	Start direction
50	0.010	2	5:00	Right
50	0.020	2	2:30	Right
50	0.040	2	1:15	Right
50	0.080	4	1:02	Right
50	0.160	5	0:37	Right
50	0.320	6	0:21	Right
- 5** New button
- 6** Update button
- 7** Remove button
- 8** ☒ Move chair to zero before test starts
- 9** ☒ Move chair to zero when test is complete
- 10** Zero pos: [°]
- Save As Default button
- Load Defaults button
- OK, Cancel, Help buttons at the bottom.

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The screenshot shows the 'Configure Sinusoidal Pendular Test' dialog box with the 'Stimulus' tab selected. The dialog has several sections for configuring the test parameters. Section 1 is 'Manual chair control' with a 'Threshold for start' of 5 [°/s]. Section 2 is 'Automatic chair control' with fields for 'Velocity [°/s]' (50), 'Frequency [Hz]' (0.32), 'Cycles' (6), and 'Acceleration: 101 [°/s²]'. Section 3 is 'Start direction left' and 'Start direction right' with a 'Test time: 0:21'. Section 4 is a table of test parameters. Section 5 is 'New', 'Update', and 'Remove' buttons. Section 6 is 'Save As Default' and 'Load Defaults' buttons. Section 7 is 'OK', 'Cancel', and 'Help' buttons.

Velocity [°/s]	Freq [Hz]	Cycles	Test time [s]	Start direction
50	0.010	2	5:00	Right
50	0.020	2	2:30	Right
50	0.040	2	1:15	Right
50	0.080	4	1:02	Right
50	0.160	5	0:37	Right
50	0.320	6	0:21	Right

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