

**Technical Note no.** : 0907  
**Date** : 2009-05-28  
**Instrument** : AS608e  
**Subject** : PC-integration/interface API

---



The AS608e version includes the Diagnostic Suite 1.0 software for transferring audiometric data to the PC for viewing, storing and printing data. For users who wish to retrieve audiometric data directly from the AS608e through USB, in order to interface/intergrate it with an external application or EMR system, we can offer a .Net-based interface/Application Programming Interface (API), including development documentation and a .Net/C# test setup project/prototype (example code). This interface requires the extended AS608e version and the Diagnostic Suite 1.0 software, i.e. the basic AS608 version does not support this feature.

Please contact Hans Pedersen ([hap@interacoustics.dk](mailto:hap@interacoustics.dk)) or technical service ([service@interacoustics.dk](mailto:service@interacoustics.dk)) for further information.