

Science **made** smarter

Why do hearing aid fitters prefer diagnostic solutions from Interacoustics?

One-stop shop, full-line supplier, service and after-sales support

Learn about products and services that will help you in your daily workflow




Interacoustics

Audiometry

Tympanometry

ABR

OAE

Hearing Aid Fitting
~~~~~

Balance



# Better fitting, happier customers

## - fewer follow ups

**A good and efficient workflow is essential when it comes to meeting the needs and expectations of your clients and making them feel confident.**

**Developed with an efficient workflow in mind**

As a leading diagnostic solutions provider, we work diligently towards improving the everyday workflow and diagnostic possibility for hearing care professionals.

This folder shows some of the products that are useful for you as a hearing aid fitter.

**The Interacoustic module concept**

Our software suites integrate various measurement modules in a common yet intuitive and user-friendly way to enable viewing, storing and printing of shared data. Protocols consisting of tests suitable for your clinic workflow can easily be set up and will optimize the fitting process.

**Service and after-sales support**

When you acquire instruments from Interacoustics, you invest in the highest quality available today. Such an investment deserves to be taken good care of and we would be pleased to assist in this regard.


Our global network of authorized distributors offer:

- Product installation
- Product training
- Calibration
- Service
- Update and installation of the latest software modules


# Integrated solutions

From video otoscopy to specialized hearing aid fitting

User-friendly integration through easy-to-use software suites



Visible speech mapping module - map the hearing instrument in real time using the Live Voice option and visually compare the response of an unaided vs. aided ear



Audiometry module - Audible Contrast Threshold (ACT™) test



Real-ear measurement module - use On-Top mode to easily adjust the hearing instruments in real time during the measurement

Workflow steps:


# Otoscopy & Tympanometry

The Viot™ video otoscope software can be integrated with the Affinity/Callisto™ Suite

The MT10 screening tympanometer for fast and easy overview of middle ear status



# Audiometry & Counseling



## Audiometry module

On top of standard air, bone and speech audiometry, the audiometry module offers a wide variety of tests.

- Audible Contrast Threshold (ACT™) test
- QuicksIN™ and TEN test (optional)
- Multiple built-in word lists

## Visible speech mapping module


Map the hearing instrument output in real time using the Live Voice option, which is very useful for counseling and managing expectations.

- Visually compare the response of an unaided vs. aided ear



# REM & HIT

Create  
your own tests  
or choose from  
the standard  
protocols



## Real-ear measurement module

Real-ear measurement made easy for a better fitting and fewer re-visits.

- Measure the output of the hearing aid in real time to ensure a proper fit
- Verification using the ISTS signal, eliminating the need to turn off the hearing instruments' adaptive features
- Use On-Top mode to easily adjust the hearing instruments in real time during real-ear measurement

## Hearing instrument testing module

The hearing instrument testing module allows for technical verification of all types of hearing instruments.

- Quickly and easily check the electro-acoustic properties of hearing aids
- Determine if the hearing aid is performing according to specifications





## Callisto™

Portable fitting solution

Callisto™ is a portable and modular PC-based system that can include both a diagnostic audiometer, a fully featured real-ear measurement module (including visible speech mapping) and a hearing instrument testing module - everything controlled from a PC.

Callisto™ accommodates the needs of the traveling clinician doing home visits or traveling between different sites.

The software concept allows you to include other Interacoustics modules within the same suite, e.g. Viot™ video otoscopy.

- Small and lightweight (565 g/1.25 lbs)
- USB powered from the PC - leave it in the bag during testing
- Advanced REM and counseling features for easy demonstration of hearing aid features
- True On-Top mode for fast and easy hearing aid fitting



TBS10 test box for hearing instrument testing with the Callisto™.



## Viot™

Video otoscope

Viot™ makes it easy to take a snapshot of your client's ear. Save it easily together with your client's data file for later use or send it to an ENT doctor for follow-up treatment.

- Ergonomic design
- Auto white balancing
- One-button control
- 90° angle view
- Windows® 10 compatible
- Full integration with Noah and OtoAccess® Database



## MT10

Fast screening tympanometer

Why initiate a full audiometric examination if a hearing loss is related to a middle ear disease? With just a few clicks, MT10 provides you with a fast status of the middle ear.

- Screening tympanometer
- Easy menu structure
- Ipsilateral reflexes (optional)
- Printout (optional)
- Data management (OtoAccess® Database and Diagnostic Suite)



## Affinity Compact

Compact and modular hearing aid fitting solution

Affinity Compact represents a new peak in hearing aid fitting - both technologically and aesthetically.

Affinity Compact combines the fast, user-friendly and well-tested Affinity Suite with a compact and modular design.

- Business-friendly modularity
- Customer-centric design
- Smart counseling
- Customizable software
- Wide range of placement options
- Hearing instrument testing
- Real-ear measurement
- Visible speech mapping
- Audiometry



Dedicated Audiometer Keyboard for the Affinity Compact.

# Science made smarter

## Interacoustics is more than state-of-the-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

[Interacoustics.com](http://Interacoustics.com)

Interacoustics A/S

Audiometer Allé 1  
5500 Middelfart  
Denmark

+45 6371 3555  
[info@interacoustics.com](mailto:info@interacoustics.com)

[interacoustics.com](http://interacoustics.com)

Go online to  
explore our  
full product  
range

## Related products



**AD528**  
Diagnostic audiometer



**AD629**  
Diagnostic audiometer



**Titan**  
Middle ear analyzer

## Product specifications

All technical and hardware specifications concerning all products can be downloaded from our website.



**Interacoustics**