

Science **made** smarter

# AD528

Diagnostics  
**made** easy

Easy to use air,  
bone & speech  
audiometer with  
superior  
ergonomics



**Interacoustics**

Audiometry  
~~~~~

Tympanometry

ABR

OAE

Hearing Aid Fitting

Balance



# Diagnostics made easy

## Air, bone & speech audiometer with superior ergonomics

### Compact & easy to use

The AD528 is the perfect audiometer for any clinic, hospital or dispensing office looking for diagnostic power that can be scaled according to your specific test needs. We proudly call it 'Diagnostics made easy':

- Starts in a few seconds. Makes it ideal for the busy clinic or traveling clinician with many power ups during a day
- Easy to use. Few buttons ensure a short learning curve
- Fast masking on/off. Just turn the right wheel to turn masking on and off

### Flexible test battery

The AD528 offers a flexible test battery. Two special tests are included. Tailor the AD528 to your needs by adding optional tests.




### Outstanding ergonomics

Ergonomics is a key aspect of the AD528. The angle of the keyboard together with the conveniently positioned tone switch and dials provides relaxed and comfortable operation for those long days of testing.

### Design your own reports and printouts

The optional software suite includes a print wizard for customizable reports and printouts. You can also print directly to a small 3" printer.



Merge and configure input from different sources in the print wizard

**Ideal for traveling**

The small footprint and low weight of 3.3lbs combined with the optional carrying bag and fast start up time make the AD528 the natural companion for the clinician on the road.

**Integrate with numerous systems**

Choose the sync license for integration to NOAH, local EMR database systems or Interacoustics OtoAccess® database.



**50 years of evolution**

Interacoustics has been developing audiometers for 50 years. This is your guarantee for a quality solution where all the benefits of previous audiometers have been compiled to match the clinical need and durability of a busy clinic that requires reliability every single day.

| Test battery                                                                                  | AD528    |
|-----------------------------------------------------------------------------------------------|----------|
| Air, bone & free-field                                                                        | x        |
| Speech (CD, MP3 , Microphone or Wave)                                                         | x        |
| Stenger                                                                                       | x        |
| Modified Hughson-Westlake                                                                     | x        |
| ABLB                                                                                          | Optional |
| SISI                                                                                          | Optional |
| Langenbeck                                                                                    | Optional |
| Békésy                                                                                        | Optional |
| Weber                                                                                         | Optional |
| Speech in noise & Binaural speech                                                             | Optional |
| QuicksIN                                                                                      | Optional |
| Sync license for integration to NOAH, local EMR systems or Interacoustics OtoAccess® database | Optional |



# Science made smarter

## Interacoustics is more than state-of-the-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

[Interacoustics-US.com](http://Interacoustics-US.com)

## Interacoustics USA

10393 West 70th Street  
Eden Prairie, MN 55344

T +1 800 947 6334  
F +1 952 903 4200

[info@interacoustics-us.com](mailto:info@interacoustics-us.com)  
[interacoustics-us.com](http://interacoustics-us.com)

Go online to  
explore our  
full product  
range

## Related products



**AT235**  
Middle ear analyzer



**AD629**  
Diagnostic audiometer



**Viot™**  
Video otoscope

## Product specifications

All technical and hardware specifications concerning all products can be downloaded from our website.



# Interacoustics