

science made smarter

Repeatable  
and hygienic  
coupler  
measurements



# Aidapter™

## Coupler measurements made consistent



**Unique Aidapter™ that replaces acoustic putty for coupler measurements.**

The Aidapter™\* from Interacoustics is a unique silicone piece that clicks onto the HA-1 adapter. You can use the Aidapter with receiver-in-canal (RIC)

and thin tube hearing aids. Compared to acoustic putty, the Aidapter™ provides a good acoustic seal and standardizes the coupling, ensuring repeatable coupler measurements. The Aidapter™ also eliminates the hassle and hygienic issues of using acoustic putty.

The Aidapter™ is a consumable good, and is compatible with all Interacoustics fitting portfolio.

\*The Aidapter™ is a patented product and all rights are owned by Interacoustics.



  
**Interacoustics**

Audiometry

Tympanometry

ABR

OAE

Hearing Aid Fitting

Balance