



# OtoAccess™ Interacoustics® Database

## User-friendly database

OtoAccess™ is a user-friendly database for storage of audiologic data and client information in private and public ENT clinics. It is equally suitable for installation on PCs as a stand alone database application as it is for networks with a number of clients using the same SQL server.

Emphasis has been placed on high-speed data handling capacity. The number of patients and sessions is virtually unlimited and very high quantities of data will have little practical impact on search speed.

OtoAccess™ has a journal feature for commenting on test sessions and a schedule planner for day-list management. Import and export facilities for XML, HL7 and GDT formats are fully implemented. A high degree of data safety ensures that OtoAccess™ meets HIPAA requirements.

Integrated HiTrack Add-In to support the collection of patient and measurement data required in newborn hearing screening programs.

## OtoAccess™ Interacoustics database

- Supports all Interacoustics products
- Rapid search & retrieval
- Unlimited storage

## Supports all Interacoustics products

### A preview of the benefits

- High speed data handling
- Unlimited patient storage
- Integrates multiple diagnostic procedures
- Journal entries for each session with standard word processing toolbar
- Import and export capabilities
- Print reports and EMR documents
- Compatible with major EMR/HIS protocols and standards, e.g. GDT, HL7 and Hi-Track
- HIPAA compliant
- Microsoft SQL Server 2005 express database for audiometric data storage
- Schedule planner with day-list
- 3 installation types: Stand-alone - Network Client - Server/Network Client



**Interacoustics**



## Interacoustics

Interacoustics is a world leading provider of diagnostic solutions in the field of hearing and balance assessment. Since 1967 we have designed and manufactured our innovative diagnostic solutions for the world of audiology with a constant focus on providing our customers with quality, dependable products.

## With you at all times

We operate in over 100 countries worldwide through a carefully selected network of distributors and service centres. Purchasing an Interacoustics product guarantees not only a quality product, but also direct access to our highly professional training and support service.

## Product specifications

All technical specifications concerning the OtoAccess™ database can be downloaded from our website.

Read more at [interacoustics-us.com](http://interacoustics-us.com)



### Interacoustics USA

10393 West 70th Street  
Eden Prairie, MN 55344

T +1 800 947 6334

F +1 952 903 4200

info@interacoustics-us.com  
interacoustics-us.com



# Interacoustics